HP Archive

This vintage Hewlett Packard document was preserved and distributed by

www. hparchive.com

Please visit us on the web !

Scanned by on-line curator: Tony Gerbic ** For FREE Distribution Only ***

4



(m) 130B





OPERATING AND SERVICING MANUAL



MODEL 130B/BR OSCILLOSCOPE

SERIALS PREFIXED: 946 -



Copyright HEWLETT-PACKARD COMPANY 1959 275 PAGE MILL ROAD, PALO ALTO, CALIFORNIA, U.S.A.

SPECIFICATIONS

Sweep Range:	0.2μ sec/cm to at least 12.5 sec/cm. 21 calibrated sweeps, accurate within ±5%, in a 1-2-5-10 sequence, 1μ sec/cm to 5 sec/cm. Vernier permits continuous adjustment of sweep time between calibrated steps and extends slowest sweep time to at least
	12.5 sec/cm.
Magnifier:	X5 Magnifier may be used on all ranges and expands fastest sweep to 0.2μ sec/cm. Accuracy within 10%.
Synchronization:	Internally from line voltage or from signals causing $1/2$ centimeter or more vertical deflection. Externally from 0.5 volts peak-to-peak or more.
Trigger Point:	Continuously adjustable from -30 to $+30$ volts on either positive or negative slope of external synchronizing signal, or from any point of the vertical signal presented on the screen.
Preset Triggering:	Switch position on sweep mode control selects optimum setting for automatic triggering.
INPUT AMPLIFIERS	Vertical and horizontal amplifiers have same characteristics.
Sensitivity:	1 mv/cm to at least 125 v/cm. 15 calibrated ranges, accurate within $\pm 5\%$, in a 1-2-5-10 sequence, 1 mv/cm to 50 v/cm. Vernier permits continuous adjustment between ranges and decreases sensitivity of 50 v/cm range to at least 125 volts/cm. Input voltage rating 600 volts dc or rms.
Phase Shift:	Within $\pm 1^{\circ}$ relative phase shift at frequencies up to 50 kc between vertical and horizontal amplifiers with verniers in cal.
Stability:	1 mv/hr after warmup.
Bandwidth:	DC Coupling: dc to 300 kc. AC Coupling: 2 cps to 300 kc. Specified bandwidth is independent of sensitivity setting.
Balanced Input:	On 1, 2, 5, 10, 20 and 50 mv/cm ranges. Cabinet Mount input impedance: 2 megohms shunted with approximately $25 \mu\mu$ f. Rack Mount input impedance: 2 megohms, approximately $125 \mu\mu$ f shunt capacity. Disconnecting the wires at the front panel which connect to the rear terminals reduces the input capacity to approximately $25 \mu\mu$ f.
Common Signal Rejection	
(Balanced input only):	Rejection at least 40 db. Common signal must not exceed 1.5 volts.

Model 130B

SPECIFICATIONS (CONT'D.)

Single Ended Input:	Cabinet Mount input impedance: 1 megohm shunted with approxi- mately 50 $\mu\mu$ f. Rack Mount input impedance: 1 megohm, approxi- mately 200 $\mu\mu$ f shunt capacity. Disconnecting the wires at the front panel connecting to the rear terminals reduces the input capacity to approximately 50 $\mu\mu$ f.
Internal Calibrator:	300 millivolts peak-to-peak $\pm 2\%$, 300 cycles squarewave applied to vertical or horizontal amplifiers by CAL position of input attenuators.
GENERAL Illuminated Graticule:	Edge lighted graticule with controlled illumination, 10 cm x 10 cm, marked in centimeter squares with 2 mm subdivisions, on major horizontal and vertical axes. Effectively shielded from ambient light.
CRT Plates:	Direct connection to deflection plates via terminals on rear. Sensi- tivity approximately 20 volts/cm.
Intensity Modulation:	Terminals on rear; 20 volts positive signal blanks CRT at normal intensity.
Cathode Ray Tube:	5 AQP mono-accelerator flat face type with 3000 volt accelerating potential. Available with P1, P2, P7 or P11 screen.
Dimensions:	Cabinet Mount: 9-3/4 in. wide, 15 in. high, 21-1/4 in. deep. Rack Mount: 19 in. wide, 8-3/4 in. high, 22 in. deep, 19-3/4 in. deep behind panel.
Weight:	Cabinet Mount: Net 41 lbs., shipping 54 lbs. Rack Mount: Net 47 lbs, shipping 62 lbs.
Power Supply:	115/230 volts ±10%, 50/1000 cycles, 160 watts.
Filter Supplied:	Color of filter compatible with screen phosphor. Green for P1 and P2, Amber for P7, Blue for P11.
Rack Mount:	Has rear terminals in parallel with front panel connections.
Accessories Furnished:	Supplied with Rack Mount. 130B-12P and Q Mounting Brackets (pr.) Two 125-57 Plugs (mate with rear terminals). Two 125-59 Clamps for 125-57 Plugs.
Accessories Available:	AC-83A Viewing Hood; face-fitting molded rubber. Additional Mounting Brackets. 130B-12P (left) and 130B-12Q (right)

Model 130B

CONTENTS

SECTION I	GENERA	AL DESCRIPTION	PAGE
	1 - 1	General	I – 1
	1 - 2	Damage in Transit	I – 1
	1 - 3	Power Line Voltages	I – 1
	1 - 4	Power Cord	I - 1
	1 - 5	Installation of Rack Mount	I - 1
		A support of the second se	
SECTION II	OPERAT	ING INSTRUCTIONS	
	2 - 1	Controls and Terminals	II - 1
	2 - 2	Rear-Access Terminals	II - 1
	2 - 3	Warm-Up Drift	II - 1
	2 - 4	AC or DC Coupling	II - 2
	2 - 5	Balanced Inputs	II - 2
	2 - 6	Operating Procedures	II - 2
SECTION III	THEORY	OF OPERATION	
	3 - 1	General Content	III - 1
	3 - 2	Over-All Operation	III - 1
	3 - 3	Vertical Amplification Channel	III - 1
	3 - 4	Horizontal Amplification Channel	III - 2
	3 - 5	Sweep Generator	III - 2
	3 - 6	Low Voltage Power Supply	III - 3
	3 - 7	High Voltage Power Supply	III - 4
	3 - 8	Calibrator	III - 4
SECTION IV	MAINTE	NANCE	
	4 - 1	Introduction	IV - 1
	4 - 2	Simple Check Procedure	IV - 1
	4 - 3	Isolating Troubles to Major Sections	IV - 3
	4 - 4	Removing the Cabinet	IV - 3
	4 - 5	Connecting for 230 Volt Operation	IV - 3
	4 - 6	Tube Replacement.	IV - 4
	4 - 7	Condensed Test and Adjustment Procedure	IV - 5
	4 - 8	Adjustment Procedure	IV - 9
	4 - 9	Turn On	IV - 9
	4 - 10	Power Supplies	IV - 9
	4 - 11	Replacing and Adjusting the CRT	IV - 11
	4 - 12	Checking and Adjusting the Calibrator	IV - 13
	4 - 13	Adjusting the Vertical Amplifier	IV - 13
	4 - 14	Adjusting the Horizontal Amplifier	IV - 15
	4 - 15	Phase Shift Adjust	IV - 15
	4 - 16	Adjusting Preset	IV - 17
	4 - 17	Adjusting the Sawtooth Generator and Sweep Amplifier	IV - 17
	100000000000000000000000000000000000000		

TABLE OF REPLACEABLE PARTS

SECTION V

- Table of Replaceable Parts . . $$\rm iv$$ 5 - 2 V - 2



CATHODE RAY TUBE WARRANTY

The cathode ray tube supplied in your Hewlett-Packard Oscilloscope and replacement cathode ray tubes purchased from ϕ , are guaranteed against electrical failure for one year from the date of sale by the Hewlett-Packard Company. Broken tubes or tubes with burned phosphor are not included in this guarantee.

Your local Hewlett-Packard representative maintains a stock of replacement tubes and will be glad to process your warranty claim for you. Please consult him.

Whenever a tube is returned for a warranty claim, the reverse side of this sheet must be filled out in full and returned with the tube. Follow shipping instructions carefully to insure safe arrival, since no credit can be allowed on broken tubes.

SHIPPING INSTRUCTIONS

- Carefully wrap the tube in 1/4" thick cotton batting or other soft padding material.
- 2) Wrap the above in heavy kraft paper.

I

I

1

I

1

1

1

1

1

ļ

1

1

1

- 3) Pack in a rigid container which is at least 4 inches larger than the tube in each dimension.
- 4) Surround the tube with at least four inches of packed excelsior or similar shock absorbing material. Be certain that the packing is tight all around the tube.
- 5) Tubes returned from outside the continental United States should be packed in a wooden box.
- Ship prepaid preferably by AIR FREIGHT or RAILWAY EXPRESS. We do not recommend parcel post or air parcel post shipment.

HEWLETT-PACKARD CO. 275 PAGE MILL ROAD, PALO ALTO, CALIF. U.S.A.

CRT WARRANTY CLAIM

() = 1×1

FROM:	DATE:
NAME:	
COMPANY:	
ADDRESS:	
Person to contact for further information	
TITLE:	
COMPANY:	
ADDRESS:	
and all the second second	
To process your claim quickly please ent	er the information indicated below:
) @INSTRUMENT MODEL	SERIAL
) TUBE TYPE	SERIAL
) ORIGINAL TUBE	REPLACEMENT TUBE
) YOUR PURCHASE ORDER NO.	and a state of the second and the second
) DATE PURCHASED	
) PURCHASED FROM	HIGH THE PARTY AND A COLOR
) COMPLAINT: (Please describe nature	e of trouble)
the man in	
) OPERATING CONDITIONS: (Please de	escribe conditions prior to and at time of failure
	SIGNATURE

SECTION I GENERAL DESCRIPTION

1-1 GENERAL

The Model 130B dc to 300 kc Oscilloscope is a general purpose oscilloscope. It can be used with either internal or external sweeps which can be either internally or externally synchronized and it can be obtained in either the cabinet or rack type mounting. Because of its high sensitivity and balanced input, the Model 130B may often be used directly with transducers, enabling you to see a direct presentation of phenomena desired without having to resort to preamplifiers.

Some of the special features of this oscilloscope are as follows:

A. LINEAR INTEGRATOR SWEEP GENERATOR

The accurate direct reading sweeps are obtained from a feedback type integrator which insures a high order linearity and stability. This type of sweep generator, is more reliable and independent of tube characteristics than other types of sweep generator.

B. X5 SWEEP EXPANSION

You speed observation and analysis of transients by expanding a two centimeter segment of the trace to 10 centimeters for easy viewing of detail. This X5 sweep expander, may be used on all sweep time settings and expands the fastest sweep time to .2 microsecond/cm.

C. CALIBRATED AMPLIFIERS

Voltage measurements of various waveforms are quickly made with the 130B, accurate within $\pm 5\%$. A built-in calibrator which is accurate within $\pm 2\%$ permits quick verification and standardization of the amplifier gain.

Phase shift measurements can be made accurately with this oscilloscope over a wide range of input frequencies.

1-2 DAMAGE IN TRANSIT

This instrument should be thoroughly inspected when it is received. If any damage is evident, refer to the "Claim for Damage in Shipment" paragraph on the Warranty sheet in this manual.

1-3 POWER LINE VOLTAGES

The Oscilloscope is shipped from the factory wired for 115 volts ac line operation, unless otherwise specified. However, the instrument may also be operated from a 230 volts ac line source if the proper conversion is made to the power transformer. This conversion is described in the Maintenance Section (Section IV).

1-4 POWER CORD

The three conductor power cable supplied with the instrument is terminated in a polarized three prong male connector recommended by the National Electrical Manufacturers' Association. The third contact is an offset round pin, added to a standard two-blade ac plug, which grounds the instrument chassis when used with the appropriate receptacle. An adapter should be used to connect the NEMA plug to a standard two contact output. When the adapter is used, the ground connection becomes a short lead from the adapter which should be connected to a suitable ground for the protection of operating personnel.

1-5 INSTALLATION OF RACK MOUNT

The 0 130BR is designed so that it can be supported in a 19 inch rack by the front panel in the usual manner; or, the dust cover may be rigidly mounted in the rack with brackets as shown in Figure 1-1. In the latter case, the chassis is supported by the dust cover and may be slipped in

Sect. I Page 2

ģ

or out easily; the screws through the front panel merely holding the chassis in place. To rack mount the 130BR using the brackets:

1) Mount the bracket as shown in Figure 1-1 with screws through the outside holes of the brackets. The length of these screws may be chosen to space the front panel from the panel rails as desired. The brackets at the rear are not necessary in most installations but can be used if added support is required. These brackets are available from the Hewlett-Packard Company as an accessory item.

2) Remove the dust cover from the 130BR and mount it in the brackets with the 10-32 trusshead screws provided.

3) Slip the 130BR into the dust cover and fasten in place with screws through the front panel.





Figure 1-1. Model 130BR Installation



SECTION II OPERATING INSTRUCTIONS

2-1 CONTROLS AND TERMINALS

Front panel operation controls are shown in Figure 2-1. This description of the operating controls enables you to operate the instrument if you have a basic knowledge of oscilloscope technique. Detailed operating procedures are given in the operating plates.

INTERNAL SWEEP CONTROLS

SWEEP TIME/CM -

This switch determines the speed at which the crt beam crosses the screen. HORIZ. SENSI-TIVITY switch must be in an INT. SWEEP position or internal sweeps are not generated. Associated with the SWEEP TIME/CM switch is a concentric VERNIER which provides continuous adjustment of sweep speed between steps. A X5 sweep magnifier operates on all ranges.

SYNC -

This three position switch lets the sweep be triggered either internally or externally. Internal triggering can be accomplished from a line frequency signal or from an applied vertical input signal of sufficient amplitude to produce a onehalf centimeter deflection. External triggering can be produced by signals having amplitude greater than 0.5 volt, peak-to-peak.

SWEEP MODE -

As this control is rotated from the extreme clockwise position, the sweep generator will pass from an un-synchronized free-running (FREE RUN) condition through a condition where only triggered operation is possible (TRIGGERED) to a position in which sweeps will not occur. At the extreme counterclockwise position the control switches into a PRESET position. This position provides optimum triggering bias for nearly all waveforms.

TRIGGER LEVEL -

This continuous control selects the level on the sync waveform where triggering is to occur. When the TRIGGER LEVEL control is set to zero, the trigger circuits are the most sensitive.

TRIGGER SLOPE -

This two-position switch, concentric with TRIGGER LEVEL, permits triggering to occur on either the positive or negative slope of internal, external or line voltage sync signals.

Horizontal or Sync INPUT -

A set of three binding posts used for receiving external sync voltages and external generated sweeping voltages. On the rack mount model only, a 3-conductor receptacle J102, mounted at the rear of the instrument, is connected in parallel with the binding posts.

2-2 REAR-ACCESS TERMINALS

-----DANGER - HIGH VOLTAGE-----

The following terminals are accessible through the rear access plate of the instrument cabinet: Horizontal and vertical deflection plates, and a terminal for crt intensity (Z-axis) modulation. See Figures 2-8 and 2-10.

2-3 WARM-UP DRIFT

When the oscilloscope is first turned on, drift in the trace will be quite noticeable, particularly at high sensitivities, the trace drift is fastest immediately following turn-on, becoming slower as the instrument warms up. Because of this drift, fine adjustment of amplifier balance should not be attempted until the instrument is thoroughly warm. For most purposes a 5 minute warm-up will be adequate.

2-4 AC OR DC COUPLING

AC coupling permits high gain to be employed without regard for the dc level involved. In the AC position the input signal (vertical or horizontal) is coupled to the amplifier through a capacitor which removes the dc component from the input. This coupling circuit has a low frequency cut-off at 2 cps. To avoid degrading input pulses or square waves below 200 cps it is advisable to use dc coupling. WHEN USING DC COUPLING THE AVER-AGE VALUE OF THE DC DETERMINES THE POSITION OF THE SWEEP ON THE OSCILLO-IF YOU ARE UNABLE TO FIND THE SCOPE. TRACE WITH THE VERTICAL POSITION CON-TROL WHEN USING DC COUPLING, TRY AC COUPLING. When AC coupled the maximum dc that may be applied is 600 volts.

2-5 BALANCED INPUTS

The instrument will accept balanced input signals on the six most sensitive ranges. This arrangement is shown in Figure 2-7.

Driving the instrument from a balanced source can be very effective in removing the unwanted stray pickup that would otherwise obscure the desired information. To take advantage of the noise reduction that is possible with a balanced input, you must be sure that neither terminal of the source is connected to ground, and use double conductor shielded cable between the source and The input cable shield must be oscilloscope. connected to a suitable ground, either at the oscilloscope or some other point. With these precautions in the external input circuit, any stray signals (noise, hum, etc.) will be coupled equally to the two input terminals, and be cancelled by the differential amplifiers. Since the desired information is applied between the two input terminals, it will be amplified and displayed in the normal manner. Since the noise is a problem

生化结果 当是一般保险的

mainly at low level, the fact that balanced input is available only on the most sensitive ranges is generally not a serious limitation.

The common-mode signal rejection will be at least 40 db (1/100 of the input signal). When using a balanced input certain limitations must be considered. The proper operating levels must be maintained on the input amplifier: The COM-MON-MODE SIGNAL VOLTAGE MUST NOT EX-CEED 1.5 VOLTS EITHER POSITIVE OR NEGA-TIVE, ON EITHER INPUT TERMINAL. Note that this is the sum of all voltages (dc plus peak ac).

NOTE

If balanced ac coupling is desired, it is necessary to connect a capacitor in the external signal path to the middle terminal, since a dc voltage on this terminal only unbalances the amplifier. This arrangement is shown in Figure 2-7.

2-6 OPERATING PROCEDURES

Basic operating procedures are described in the following illustrations. Positions of controls are different on the cabinet model but their functions are identical to those of the rack model.

Figure	Description		

- 2-2 VERTICAL BALANCE ADJUSTMENT
- 2-3 HORIZONTAL BALANCE ADJUSTMENT
- 2-4 INTERNAL SWEEP-INTERNAL SYNCHRO-NIZATION
- 2-5 INTERNAL SWEEP-EXTERNAL SYNCHRO-NIZATION
- 2-6 EXTERNAL HORIZONTAL INPUT
- 2-7 AC COUPLING BALANCED INPUT
- 2-8 CONNECTION TO CRT DEFLECTION PLATES
- 2-9 EXTERNAL INTENSITY MODULATION
- 2-10 ALIGNING SCOPE TRACE WITH GRATICULE



- 4. Short vertical input terminals together.
- 5. Set AC-DC switch to DC.

factory.)

- 6. Turn VERNIER to CAL.
- 7. Center trace with VERT. POS. control.

- Set VERT. SENSITIVITY to 5 MILLIVOLTS/ CM, repeat step 9.
- 11. Set VERT. SENSITIVITY to 1 MILLIVOLT/ CM, repeat step 9.
- 12. Turn VERNIER fully CCW, repeat step 9.

Repeat steps 6 through 12 until the trace does not move as the VERNIER control is rotated.

In Figure 2-3, delete step 9. A separate adjustment (Bal. Adj. control on the etched board) is provided to balance the



JUMPER

1

After Warm-Up:

- 1. Short together the horizontal input terminals.
- 2. Set AC-DC switch to DC.
- 3. Set HORIZ. SENSITIVITY to 50 MILLI-VOLTS/CM.
- 4. Turn VERNIER to CAL.
- 5. Center spot with HORIZ. POS. control, if necessary, adjust vertical position.
- 6. Set HORIZ. SENSITIVITY to 5 MILLI-VOLTS/CM.

- 7. Return the spot to center with coarse (screwdriver) HORIZ. DC BAL. control or with the fine (knob) control if the unbalance is slight.
- 8. Set HORIZ, SENSITIVITY to 1 MILLIVOLT/ CM, repeat step 7.

Repeat steps 4 through 9 until the spot does not move as the VERNIER control is rotated.

Figure 2-3

. 130B

Model 130B



- Set HORIZ. SENSITIVITY switch to INT. SWEEP X1 (or to X5 for magnified sweeps)
- 2. Set SYNC switch to INT.
- 3. Set SWEEP MODE to PRESET.
- 4. Connect vertical input signal into vertical input terminals.
- 5. Set AC-DC switch for type coupling desired.
- 6. Adjust VERT. SENSITIVITY for desired sensitivity.

- Set TRIGGER SLOPE switch for triggering on positive or negative slope of input signal, as desired.
- 8. Set TRIGGER LEVEL control to 0.
- 9. Select desired sweep speed with SWEEP TIME/CM switch.
- 10. Adjust TRIGGER LEVEL to start trace at desired level. In some cases, it may be necessary to switch SWEEP MODE from PRESET to an individual adjustment for the particular trace being viewed.

Sect. II Page 6

Model 130B



- Set HORIZ. SENSITIVITY switch to INT. SWEEP X1 (or to X5 for magnified sweeps).
- 2. Set SYNC switch to EXT.
- 3. Set SWEEP MODE to PRESET.
- 4. Feed synchronizing signal (0.5 volts p-p or more) to the horizontal input terminals.
- 5. Set AC-DC switch for type coupling desired.
- 6. Set TRIGGER LEVEL to 0.
- 7. Feed vertical input signal into vertical input terminals.

- 8. Adjust VERT. SENSITIVITY for desired sensitivity.
- 9. Select desired sweep speed with SWEEP TIME/CM switch.
- 10. Set TRIGGER SLOPE for triggering on positive or negative slope, as desired.
- 11. Adjust TRIGGER LEVEL to start trace at desired level. In some cases, it may be found necessary to switch SWEEP MODE from PRESET to an individual adjustment for the particular trace being viewed.

EXTERNAL HORIZONTAL INPUT



- 1. Feed horizontal signal to horizontal input terminals.
- 2. Set AC-DC switch for type of input coupling desired.
- 3. Set HORIZ. SENSITIVITY switch for desired sensitivity.
- 4. Adjust horizontal position of pattern with HORIZ. POS. control.

This type of input will be found useful for viewing Lissajous patterns, etc.



The following procedure is for the vertical input, but is the same for the horizontal input.

- Set VERT. SENSITIVITY to 50 MILLIVOLTS/ CM, input not balanced on higher ranges.
- 2. Set AC-DC switch to AC.
- 3. Disconnect shorting strap.

- 4. Connect 0.1 microfarad capacitor to midterminal.
- 5. Connect input signal to A and B.
- 6. Ground input at the black terminal.

The capacitor must be used to block any dc.



The following procedure is for connecting external signals to the vertical deflection plates, but is the same for the horizontal plates.

- 1. Remove rear access plate fastened by four screws.
- 2. Remove the shorting bars between the Vertical Amplifier and terminals D3 and D4 and replace them with 1 megohm, 1/2 watt resistor.

For balanced AC coupling:

3. Connect balanced signal through appropriate capacitor to D3 and D4.

For single-ended AC coupling:

- 4. Bypass D4 to chassis with an adequate capacity.
- Connect the signal to D3 through an appropriate capacitor.

NOTE: If it is desired to have positive voltage deflect the beam downward, bypass D3 to chassis and connect the signal to D4.

Sect. II Page 10





To intensity modulate the CRT with external signals:

- 1. Remove rear access plate fastened by four small screws at rear of dust cover.
- 2. Remove shorting bar.
- 3. Connect modulating signal to these terminals. A positive voltage of 20 volts peak will blank the CRT trace from normal intensity.

ALIGNING SCOPE TRACE WITH GRATICULE



CAUTION: DANGEROUS VOLTAGES ARE PRESENT INSIDE THE INSTRUMENT

Remove two screws at rear of dust cover and slide cover off to rear. Fiber lever (2) controls both radial and longitudinal positioning of CRT and is locked by clamp (1).

To align sweep trace with graticule loosen clamp (1) with a screwdriver. Rotate fiber arm (2) until the trace is parallel to horizontal lines on graticule. Tighten clamp (1) after adjustment has been made.

AC VERT. SENSITIVITY QDC 0 VERTICAL INPUT 0 INPUT INTER-STAGE INPUT OUTPUT ATTENUATOR ATTENUATOR AMPLIFIER AMPLIFIER -> \rightarrow 0 0 (SINGLE-ENDED) (BALANCED) CAL. ← 6.3 VAC FROM T30 O SWEEP TIME / CM LINE Ó TRIGGER SAWTOOTH CALIBRATOR INT. GENERATOR GENERATOR SYNC. EXT. AC SYNC OR HORIZONTAL QDC INPUT 0 0 C INPUT INTER-STAGE INPUT OUTPUT ATTENUATOR ATTENUATOR 0 AMPLIFIER AMPLIFIER (SINGLE-ENDED) (BALANCED) RO 0 0 CAL. HORIZ. SENSITIVITY 0

Figure 3-1. Model 130B Block Diagram

Model 130B

Sect. III Page 0

SECTION III THEORY OF OPERATION

3-1 GENERAL CONTENT

This section contains a brief description of the over-all operation of the Model 130B Oscilloscope, description of each major section and detailed description of a Schmitt trigger.

3-2 OVER-ALL OPERATION

The block diagram in Figure 3-1 shows the basic circuits of the Model 130B Oscilloscope.

A. VERTICAL AMPLIFIER

The Vertical Amplifier receives the input signal, amplifies it, and drives the vertical deflection plates of the cathode ray tube. In addition, this amplifier determines the vertical position of the spot on the screen and supplies a signal for synchronizing the sweep with the vertical input signal.

B. HORIZONTAL AMPLIFIER

The Horizontal Amplifier receives its signal either from the horizontal INPUT jack or from the Sweep Generator, amplifies it and drives the horizontal deflection plates of the cathode ray tube. Except for the provisions in the Horizontal Amplifier for amplifying the internally-generated sawtooth voltage, it is essentially the same as the Vertical Amplifier.

C. SWEEP GENERATOR

The Sweep Generator forms a sawtooth voltage to control the horizontal movement of the spot across the face of the cathode ray tube. The Sweep Generator is divided into two parts: 1) a sawtooth generator, 2) a trigger generator, which starts the sawtooth. The trigger generator controls allow the operator to choose the point at which the sawtooth sweep begins. In addition to forming the internal sweep of the oscilloscope, the Sweep Generator also supplies the required unblanking pulse which brightens the trace during each sweep.

D. CALIBRATOR

An internal square-wave calibrator, with a nominal frequency of 300 cps, is provided for setting the basic gain of the amplifiers. Turning either the VERT. or HORIZ. SENSITIVITY switches to CAL., turns on the calibrator supply voltage and connects its output to the appropriate amplifier.

E. CATHODE RAY TUBE

The cathode ray tube is a 5AQP - monoaccelerator type. It is normally supplied with a P1 phosphor screen but is available in the P7 and P11 phosphors also and P2 upon special order. All are electrically interchangeable and the tube is easily changed. The mono-accelerator anode makes possible a simple astigmatism adjustment which requires no resetting when adjusting the FOCUS or INTENSITY controls. The deflection plate terminals are connected through removable jumpers at the rear of the instrument so that direct connections to the plates can be made easily.

3-3 VERTICAL AMPLIFICATION CHANNEL

The vertical amplification channel consists of three parts: the AC-DC switch, the input attenuator, and the amplifier section proper.

A. AC - DC SWITCH

The signal comes into the input terminals and is fed to the AC-DC switch. For ac coupling, a capacitor is switched into the signal path. In the DC position, the signal goes directly to the input attenuator.

B. INPUT ATTENUATOR

The input attenuator is a sixteen position switch having fifteen calibrated ranges (1 MILLI-VOLT/CM to 50 VOLTS/CM) and a calibrate position. When the switch is in the CAL, position, the input of the amplifier is directly connected to the output of the internal calibrator. On ranges less sensitive than 50 MILLIVOLTS/CM, singleended frequency-compensated attenuators are inserted ahead of the Vertical Amplifier. On the six most sensitive ranges, balanced-type attenuators are inserted between the second differential amplifier (V2) and the third differential amplifier (V3). On the six most sensitive ranges, balanced input signals may be applied to the input terminals after removing the jumper to the ground terminal. The sensitivity may be varied continuously between ranges by means of the VERNIER control.

C. VERTICAL AMPLIFIER

The Vertical Amplifier consists of four stages of balanced differential amplifiers* in cascade. The first stage (V1) has the VERT, DC BAL, adjustment (R10A, B) in its cathode circuit which adjusts the current division between the two halves of the stage. The second stage has a VERNIER control in the cathode circuit which varies the gain of the amplifier between ranges of the VERT. SENSITIVITY switch, and another dc balance adjustment (R20) is also provided. In the last three stages, neutralizing capacitors are used to cancel the coupling effects between the input and output of the amplifier arising from the inter-electrode capacitances. The output of the second stage is fed to the balanced attenuator of the VERT, SEN-SITIVITY switch. The output of the balanced attenuator is connected to the third balanced differential amplifier (V3). The third stage has two potentiometers in its cathode circuit, one controls the vertical position of the pattern (VERT, POS) and the other adjusts the basic gain of the Vertical Amplifier (R40, Gain Adj.). The fourth balanced differential amplifier (V4) is the output stage. The neon lamps in the grid-cathode circuit of V4 protect the tube when the Model 130B is first turned on. The output of V4 drives the vertical deflection plates of the cathode ray tube. In addition, synchronization signals are coupled from

* Valley and Wallman, "Vacuum Tube Amplifier", Massachusetts Institute of Technology Radiation Series, vol. 18, pp 441-451. McGraw-Hill Book Company, Inc., New York, 1948. the plates of V4 and coupled into the Sweep Generator to trigger the sweep during either IN-TERNAL + or INTERNAL - synchronization. As a precaution against drift and hum, a regulated dc supply is used for the heaters of the first three stages.

3-4 HORIZONTAL AMPLIFICATION CHANNEL

The Horizontal Amplifier is essentially identical to the Vertical Amplifier, except in the INT. SWEEP X1 and X5 position of the HORIZ. SEN-SITIVITY switch. In these positions, the sawtooth signal from the Sweep Generator is fed through the sweep attenuator to the grid of V103, the third balanced differential amplifier. In the INT. X5 position, R164, X5 Mag. Adj., in the cathode circuit of V104 sets the gain of the amplifier to obtain sweep magnification of X5. The output of V104 drives the horizontal deflection plates of the cathode ray tube.

3-5 SWEEP GENERATOR

The sweep generator provides a sawtooth voltage to produce linear horizontal movement of the spot across the face of the cathode ray tube when the HORIZ. SENSITIVITY switch is set to INT. SWEEP (X1 or X5). In addition, the sweep generator furnishes the pulse required to unblank the cathode ray tube during each sweep.

The sweep generator consists of a Trigger Generator, a Sawtooth Generator, and a Gate Out Cathode Follower.

A. TRIGGER GENERATOR

The purpose of the Trigger Generator is to receive a synchronizing signal and convert it into a fast, constant-amplitude pulse to start the Sawtooth Generator.

The Trigger Generator consists of a SYNC selector switch (S201), a Trigger Amplifier (V201), and a Trigger Generator (V202). The SYNC selector switch accepts a signal from:

- the Vertical Amplifier (internal synchronization, + or -),
- 2) an internal 6.3 volt source (line-frequency synchronization), or
- the horizontal INPUT terminals (external synchronization).

The synchronizing signal is fed to V201 which amplifies the signal and delivers it in the proper phase, as selected by the TRIGGER SLOPE switch, to the Trigger Generator. Adjustment of the TRIGGER LEVEL control sets the output level of V201, determining the point on the input waveform that will trigger the Trigger Generator (V202). Trigger Generator (V202) is a Schmitt trigger circuit; a discussion of the Schmitt trigger follows:

A Schmitt trigger consists of two amplifiers, A and B, having both plate-to-grid and cathode-to-The circuit has two stable cathode coupling. states: A side conducting, B side cut off; B side conducting, A side cut off. Due to regenerative action the change-over from one state to the other is very rapid, producing fast rise and decay times in the square-wave output. The levels at which the change-over takes place (hysteresis limits) can be adjusted to be close together as in the Trigger Generator (V202) or widely spaced as in the Start-Stop Trigger (V203). To trigger the circuit, the A side grid voltage must cross a particular hysteresis limit to change the state of the circuit. For example, if the A side is conducting, driving the grid voltage positive through the upper hysteresis limit will have no effect, but driving the grid voltage negative through the lower hysteresis limit will put the A side out of conduction and B side into conduction.

B. SAWTOOTH GENERATOR

The Sawtooth Generator consists of Start-Stop Trigger (V203), and Integrator Switch (V205), a Feedback Integrator (V206B), and Integrator Cathode Follower (V206A), and a Retriggering Hold-Off Cathode Follower (V207B).

Start-Stop Trigger (V203), a Schmitt trigger circuit, is fed by Trigger Generator (V202). The square wave output of V203 is fed directly to the Integrator Switch (V205), which in turn controls the action of Feedback Integrator (V206B). When V203 produces a negative pulse, it causes V205 to cut off permitting V206B to commence operation.

Feedback Integrator (V206B), a Miller integrator circuit^{*}, generates essentially a positive linearly rising waveform, which is applied to the Horizontal Amplifier to sweep the trace across the face of the cathode ray tube. The rate at which this sweep takes place is determined by the values of the RC network in the grid circuit of V206B. These values are varied by the SWEEP TIME switch. The output of V206B is fed through a neon lamp (I203) to the Integrator Cathode Follower (V206A). Neon lamps (I204 to I206) are used to drop the voltage to the proper level and at the same time furnish a direct-coupled path for the signal. The neon lamps are shunted with a capacitor to improve the high-frequency response of the circuit, and a series resistor is used to eliminate any tendency toward oscillation.

The output of the Integrator Cathode Follower (V206A) is fed to two circuits: 1) through the sweep attenuator to the Horizontal Amplifier and 2) to the Retriggering Hold-Off Cathode Follower (V207B) in the Sawtooth Generator feedback circuit. During the Sweep, V207B conducts and the capacitor in its cathode circuit charges. However, at the termination of the sweep, V207B is cut off and the cathode capacitor discharges, maintaining a positive bias on the grid of V203A. This hold-off bias allows sufficient time between sweeps for the Sweep Generator to recover. The bias which determines the triggering level of the Start-Stop Trigger (V203A) is supplied by the Retriggering Bias Control (V207A). The bias is adjusted by the SWEEP MODE control, R218, in the grid circuit of V207A.

C. GATE OUT CATHODE FOLLOWER

Another function of the Start-Stop Trigger is to furnish a pulse to unblank the cathode ray tube. The Gate Out Cathode Follower (V204), couples the required positive unblanking pulse from the Start-Stop Trigger to the grid of the crt for the duration of the sweep.

3-6 LOW VOLTAGE POWER SUPPLY

The low-voltage power supply consists of four regulated voltage supplies, three positive (+585V, +300V, +100V) and one negative (-150V), furnishing the plate voltages and dc filament voltages required for the instrument.

The operation of each of the four regulators is similar; only the -150 volt supply will be discussed. V306, V307 and V308 constitute the voltage regulator circuit for the -150 volt supply. V308, a glow discharge tube, provides a fixed reference voltage for the cathode of V307, the

^{*} Millman and Taub, "Pulse and Digital Circuits" pp 216-228, McGraw-Hill Book Company, Inc., New York, 1956.

Sect. III Page 4

Control Tube. V306 operates as the regulator tube (or variable resistor), controlled by the voltage at the plate of V307. If the regulated output from the cathode of V306 tends to increase, the voltage at the grid of V307 tends to increase, causing V307 to draw more current. This lowers the plate voltage of V307 and the grid voltage of V306, resulting in greater plate resistance for V306. Increased plate resistance causes a greater voltage drop across V306, compensating for the increased output voltage from the regulator and resulting in a substantially constant output.

If the regulated output tends to decrease, the reverse of the above action occurs, tending to maintain the output voltage constant. In the same manner, changes in voltage at the screen of V307 compensate for variations and ripple in the input (unregulated) voltage. The output of the -150 volt supply serves as the reference voltage for the three positive-voltage supplies.

3-7 HIGH-VOLTAGE POWER SUPPLY

The high-voltage power supply provides regulated dc voltage to the cathode and control grid of the cathode ray tube. The high-voltage power supply consists of an RF Oscillator tube (V313), a highvoltage transformer (T302), high-voltage rectifiers (V310,311) and a High-Voltage Control Tube (V312). The RF Oscillator, a Hartley circuit, oscillates at a frequency of approximately 100 kc. The high-voltage transformer has two separate secondaries which feed the High-Voltage Rectifiers.

The output of V310 is connected to the cathode of the cathode ray tube. A fraction of this voltage is fed to the High-Voltage Control Tube V312, a dc-coupled amplifier. The output of V312 is fed back to the screen of RF Oscillator tube (V313) in proper phase to oppose any change in the highvoltage output. The INTENSITY control in the output of this supply determines the voltage on the cathode of the cathode ray tube.

The output of V311 is connected to the control grid of the cathode ray tube, and normally the crt beam is cut off. During the sweep operation, a positive pulse from the Gate Out Cathode Follower (V204) in the Sweep Generator circuit overrides the negative crt grid cutoff voltage and unblanks the cathode ray tube. The brilliance of the trace may be adjusted with the Intensity Adjust potentiometer (R343), in series with grid-voltage supply.

3-8 CALIBRATOR

The Calibrator, a square-wave oscillator, produces an accurate voltage across R244 for application to either amplifier for setting the basic gain. Turning either the VERT. or HORIZ. SENSITIVITY switches to CAL. turns on the Calibrator and connects its output to the appropriate amplifier.

The Calibrator consists of two neon lamps (1207 and 1208) in a relaxation oscillator circuit. Operation of the Calibrator is as follows:

When the +300 volt supply is applied to the Calibrator, 1207 will ionize first due to higher potential across it compared to the potential across 1208. When 1207 fires it will draw current through R243. However, the voltage at the junction of R242, C213 and R243 will build up slowly because the voltage across a capacitor cannot change instantaneously. As C213 allows this voltage to change, the voltage at the common junction of I 207 and I 208 will also change, since the voltage drop across the ionized neon lamp is constant (approximately 60 volts). As the voltage at the common junction of I 207 and I 208 reaches approximately +70 volts, 1208 will fire. This additional current through R240 and R241 will reduce the voltage across I 207 and it will de-ionize. I 208 remains lit until the voltage across C213 charges through R243 to a voltage approximately 70 volts below the voltage that appears at the common junction of I 207 and I 208. I 207 will now fire and the action will repeat itself.

I 208 is thus alternately turned off and on at a rate of approximately 300 cps. The output of the Calibrator is taken from the current passing through R244 and I 208. The output is approximately a square wave which can be set with R240 to obtain 300 millivolts in amplitude.

SECTION IV

4-1 INTRODUCTION

This section contains instructions for testing, adjusting, and trouble shooting the Model 130B Oscilloscope.

Standard, readily available components are used for manufacture of @ instruments whenever possible. Special components are available through your local @ Representative who maintains a part stock for your convenience.

When ordering parts, specify instrument model and serial number plus the component description and stock number appearing in the Table of Replaceable Parts.

Your local Φ Representative maintains complete facilities and specially trained personnel to assist you with any problems you may have with Φ instruments.

The material in this section is divided according to circuit functions, each section having a complete set of adjustment instructions. The material in this section is as follows:

- 4-2 Simple Check Procedure
- 4-3 Isolating Troubles to Major Sections
- 4-4 Removing the Cabinet
- 4-5 Connecting for 230 Volt Operation
- 4-6 Tube Replacement
- 4-7 Condensed Test and Adjustment Procedure
- 4-8 Adjustment Procedure
- 4-9 Turn On
- **4-10** Power Supplies
- 4-11 Replacing and Adjusting the CRT
- 4-12 Checking and Adjusting the Calibrator
- 4-13 Adjusting the Vertical Amplifier
- 4-14 Adjusting the Horizontal Amplifier
- 4-15 Phase Shift Adjust
- 4-16 Adjusting Preset
- 4-17 Adjusting the Sawtooth Generator and Sweep Amplifier

The following test equipment is used for testing and adjusting the Model 130B Oscilloscope during manufacture. Equivalent test equipment may be used.

1) A high impedance dc vacuum tube voltmeter, such as an $\frac{1}{9}$ Model 410B with an $\frac{1}{9}$ Model 459A DC Voltage Multiplier, calibrated to an accuracy of $\pm 1\%$.

2) A high impedance ac vacuum tube voltmeter, such as an $\frac{1}{20}$ Model 400D/H/L.

3) A variable power line transformer with a minimum rating of 3 amps.

4) A square-wave generator such as an @ Model 211A.

5) A sine-wave oscillator with a maximum frequency of at least 500,000 cycles, such as an $\frac{1}{20}$ Model 200CD.

6) An accurate time mark generator suitable for sweep speed calibration.

4-2 SIMPLE CHECK PROCEDURE

This check should be performed first whenever instrument malfunction is suspected. It is not necessary to remove the instrument from the cabinet.

Set both VERT, and HORIZ. SENSITIVITY switches on CAL. The pattern should be a straight line tilted at 45 degrees. In addition, the deflection should be a total of six centimeters in the horizontal and vertical directions.

If the proper pattern is obtained, it is likely that both the Vertical and Horizontal Amplifier, the Power Supplies and the Calibrator are functioning properly. To check the Sweep Generator proceed as follows: Sect. IV Page 2



1) Leaving the VERT. SENSITIVITY switch in CAL, switch HORIZ. SENSITIVITY switch to INT. SWEEP X1.

2) Switch SWEEP TIME/CM switch to 1 MILLI-SECOND. A six centimeter square-wave pattern should appear on the screen. If no pattern is obtained be sure SWEEP MODE is in PRESET, SYNC switch is in INT., and adjust TRIGGER LEVEL to trigger. If a pattern cannot be obtained, the malfunction is most likely in the Sweep Generator.

4-3 ISOLATING TROUBLES TO MAJOR SECTIONS

Determining which major section contains a malfunction is usually not a difficult process, if the following general rules are remembered.

1) A failure affecting all major sections can usually be traced to the power supply.

2) A failure occurring in the last two stages of the Horizontal Amplifier also will affect internally generated sweeps, while a failure in the first two stages affects only the Horizontal Amplifier.

3) A sweep Generator failure affects internally generated sweeps only, and does not affect the Horizontal Amplifier.

4) If following the Simple Check Procedure does not produce a trace or spot on the screen, measure the voltages on the deflection plates of the Cathode-Ray Tube (deflection plate terminal board is a convenient place to measure). If, with both VERT, and HORIZ, SENSITIVITY switches set to 50 MILLIVOLTS/CM, these voltages can be set to approximately 480 vdc using the position controls, look for trouble in the high voltage section of the power supplies. If one set of deflection plates has unbalanced voltages, or if the position control must be turned far from its mechanical center to balance these voltages, look for trouble in that amplifier. If both sets of deflection plates have unusual voltages, look for trouble in the power supply.

5) If the series heater string should open, all major sections will be inoperative.

6) The two sides of the direct-coupled differential amplifier, such as are used in the Vertical and Horizontal Amplifiers on the 130B, are balanced and, unless a signal is present, the spot will be motionless in the center of the screen. Any signal,

whether this signal is applied to the input terminals or is supplied by an internal source, such as a positioning or balance control, causes the spot to move from the center of the screen. As the instrument ages it is to be expected that a drift will occur which must be compensated by internal adjustments. However, should there be a component failure in either amplifier the spot will be thrown off the screen and usually out of range of adjustment of the balance and positioning controls. To isolate the trouble, begin by shorting together the grids of the amplifier closest to the output. If the trace (spot) returns to the screen, the fault is ahead of this stage. Proceed towards the front of the amplifier. If shorting the grids of one stage does not return the spot to the screen, the fault is in this stage, or if there is a balancing control in this stage, it may be out of adjustment.

7) To check the Sweep Generator quickly, set the SWEEP TIME/CM switch to 5 or 10 MILLI-SECONDS/CM, turn the SWEEP MODE control to FREE RUN, and observe 1201, 1202 and 1203. These are the three neon lamps near V206 (6AW8) on the Sweep Generator etched circuit board. If these lamps flicker regularly, the Sweep Generator is sweeping. Turning the SWEEP MODE into the TRIGGER region should stop the generation of sweeps and, hence, the flickering of the neon lamps.

4-4 REMOVING THE CABINET

In the cabinet model, remove the two screws at the rear of the cabinet, and push the instrument forward.

If the 130BR has been rack-mounted with brackets as described in Figure 1-1, remove the screws which pass through the front panel, and withdraw the chassis. If the instrument is out of the rack, turn it on its face (handles will protect the controls), remove the two screws at the rear, and lift off the dust cover.

4-5 CONNECTING FOR 230 VOLT OPERATION

Unless otherwise requested by the customer, Φ instruments are shipped with their power transformer primaries connected in parallel for operation on 115 volt (nominal) power lines.

To convert to 230 volt supply, remove the instrument from its cabinet or dust cover by removing the two screws at the rear of the chassis, and Sect. IV Page 4

push the chassis forward. At the primary of the power transformer (marked A), remove the wires connecting terminals 2 and 5, and 1 and 4. Then connect 1 to 2 as shown in Figure 4-2, and replace the 2 amp slow-blow fuse (F301) with a 1-1/4 amp slow-blow fuse. The instrument may now be connected to the 230 volt line.



Figure 4-2. Line Voltage Connection

BAUADO MELY I WHE RECEIPT HAR STITLE A TOMA

4-6 TUBE REPLACEMENT

In many cases instrument malfunction can be corrected by replacing a weak or defective tube. Before changing the setting of any internal adjustment, check the tubes. Adjustments made in an attempt to compensate for a defective tube will often complicate the repair problem.

It is a good practice to check tubes by substitution rather than by using a "tube checker". The results obtained from the "tube checker" can be misleading. Before removing a tube, mark it so that if the tube is good it can be returned to the same socket. Replace only tubes proved to be weak or defective.

Any tube with corresponding standard ElA (JEDEC) characteristics can be used as a replacement. Where variation in tube characteristics will affect circuit performance, an adjustment is provided. The following table lists the tests and adjustments which should be performed if such tubes are replaced.

The chart in Table 4-2 lists all tubes in the 130B with their functions and adjustments required when replacing tubes. The heaters of some tubes are operated in series from a regulated dc voltage obtained from the Low-Voltage Power Supply. These tubes are identified in the chart with an asterisk and their heaters are shown in the Filament and Primary Detail Schematic. If a tube in the dc string is pulled or burned out, all tubes in the string will be turned off.

4-7 CONDENSED TEST AND ADJUSTMENT PROCEDURE

All basic tests and adjustments are covered in the following Table 4-1. In most cases, this table will cover all normal adjustment needs for the oscilloscope. For a more complete and detailed test procedure refer to paragraph 4-8.

If the instrument is not operating, refer to paragraphs 4-3 and 4-6.

If a tube is replaced, refer to Table 4-2 and complete the indicated adjustments.

		And the second second second second second					
TTATT T	4 1	CONDENGED	TOTOOT	ANT	A TO TELEVIER A LEADER	DDOOEDIDE	
TABLE	4-1.	CONDENSED	ILOI	AND	ADJUSIMENT	PROCEDURE	
			100 Mar 1997		Strengthered and the strengthere and the s		

	Test	External Equip- ment Required	Procedure	Adjust	Notes
1.	Low Voltage Power Supply	DC vtvm with 1% accuracy	Measure all low voltage power supply outputs should be within the following limits: -150 \pm 6 volts $+100$ \pm 4 volts $+300$ \pm 12 volts \pm 585 \pm 25 volts	If voltages are outside limits, adjust R332 for -150 volts.	Check sweep calibration if -150V is adjusted.
2.	Vertical ampli- fier balance	NONE	HOR.SENS. to INT.SWEEP X1, SWEEP MODE to free-run, SYNC to INT., SWEEP TIME to 1.0 ms/cm, short-circuit input terminals and set INPUT for DC. VERT.SENS. and Vernier to Cal.	Center bottom of calibrating signal with VERT.POS. control	Repeat as required.
			VERT.SENS. to 1 mv/cm, Center VERT.DC.BAL control (knob).	Center trace with coarse bal- ance control. (Screw adjustment in center of DC BAL control.)	
3.	Vertical VERNIER balance	NONE	Short circuit input terminals and set INPUT for DC. VERT.SENS. to 1 mv/cm, VERNIER to Cal.	Center spot (or trace) with VERT. POS. control.	Repeat as required.
1		and the second	VERNIER fully CCW	Return spot to center with R20	
4.	Vertical ampli- fier gain	400 cycle Voltage Cali- bration Gen- erator	VERT.SENS. to 50 mv/cm. VERNIER to Cal. Connect 300 mv p-p from Calibration Gen- erator to vertical input.	Adjust R40 for 6 cm deflection.	
		Square Wave Generator	Connect 50 kc square wave to Vert. Input. Adjust square wave generator for 6 cm de- flection. SYNC to INT, Adjust SWEEP MODE and TRIGGER LEVEL for stable picture.	Adjust C12 for best square wave.	
5.	Calibrator	NONE	VERT.SENS. and VERNIER to CAL, SWEEP MODE to free- run; SWEEP TIME to 1 ms/cm.	Adjust R240 for 6 cm deflection.	
6.	6. Horizontal amplifier balance	NONE	Short-circuit input terminals and set INPUT for DC. VERT.SENS. to 50 mv/cm, with no input. HOR.SENS. and VERNIER to CAL.	Center the left spot with the HOR.POS. control.	Repeat as required.
			HOR.SENS. to 1 mv/cm, Center the HOR.DC BAL con- trol (knob).	Center the spot with the coarse balance control (screw driver adjustment in center of DC BAL control).	
7.	Horizontal VERNIER balance	NONE	Short-circuit input terminals and set INPUT for DC. HOR.SENS. to 1 mv/cm, VERNIER to CAL.	Center spot with POS. control HORIZ.	Repeat as required.
	_		VERNIER fully CCW.	Return spot to center with R120.	
8.	Horizontal amplifier	NONE	HOR.SENS. and VERNIER to CAL.	Adjust R144 for 6 cm between spots.	
	gain	Square wave generator	HOR.SENS. to 50 mv/cm, Connect 50 kc square wave to Hor. input and adjust for 6 cm deflection.	Adjust C114 for best defined spots.	
9.	Sweep gain	Time Marker Generator	HOR.SENS. to INT.SWEEP X1 SWEEP TIME to 1 ms/cm, VERNIER to CAL; 1 kc markers from generator to VERT. input. SYNC to INT. Adjust SWEEP MODE and TRIG.LEVEL for stable pattern.	Adjust R134 for one marker/cm	
			HOR.SENS. to INT.SWEEP X5.	Adjust R164 for markers 5 cm apart.	
10.	Sweep preset	DC VTVM	HOR.SENS. to INT.SWEEP X1. SWEEP TIME to 1 ms/cm, SWEEP MODE to PRE-SET, SYNC to EXT. with no input. Connect VTVM 30 volt range between center arm of Preset pot (R220) and ground.	Slowly adjust R220 and note volt- age just prior to sweep start. Adjust pre-set for 2 volts more positive than voltage noted.	
11.	Sweep length	Sine Wave Oscillator	HOR.SENS. to INT.SWEEP X1. SWEEP TIME to 1 ms/cm. SYNC to INT. Connect 500 kc sine wave to vertical input. Adjust level and VERT SENS. to produce 6 cm vertical deflection.	Adjust R229 for a trace about 10.5 cm long.	

Sect. IV Page 7

Model 130B

TABLE 4-2. TUBE REPLACEMENT CHART

Ref.	Tube	Function	Adjustment
description of		VERTICAL AMPLIFIER	
V1* V2* V3* V4*	12AU7‡ 12AU7 12AT7 6DJ8/6BQ7	Phase Inverter Amplifier Differential Amplifier Differential Amplifier Differential Amplifier	Vertical Amplifier (par. 4-13A/B) Vertical Amplifier (par. 4-13A/B) Vertical Amplifier (par. 4-13B) Vertical Amplifier (par. 4-13B)
		HORIZONTAL AMPLIFIER	
V101* V102* V103* V104*	12AU7‡ 12AU7 12AT7 6DJ8/6BQ7	Phase Inverter Amplifier Differential Amplifier Differential Amplifier Differential Amplifier	Horizontal Amplifier (par. 4-14A/B) Horizontal Amplifier (par. 4-14A/B) Horizontal Amplifier (par. 4-14B) Sawtooth Generator (par. 4-17)
		SWEEP GENERATOR	
V201* V202 V203	6DJ8/6BQ7A 12AT7 6U8	Trigger Amplifier Trigger Generator Sweep Start-Stop Trigger	none none Adj. Preset, Sweep Length
V204 V205*	6C4 12AL5	Gate Out Cathode Follower a. Integrator Switch	(par. 4-15 & 4-10) none
V 206	6AW8	 b. Integrator Switch a. Integrator Cathode Follower b. Feedback Integrator 	none
V207	12AX7	a. Retriggering Hold Off b. Retriggering Bias Control	none
		POWER SUPPLY	
V301 V302 V303 V304*	12B4 6AU6 12B4 6BH6	+300-volt Series Regulator +300-volt Control Tube +100-volt Series Regulator +100-volt Control Tube	none none none none
V305 V306 V307* V308	6X4 12B4 6BH6 5651	-150-volt Rectifier -150-volt Series Regulator -150-volt Control Tube Reference Tube	none none LV Supply (par. 4-10A)
V309 V310 V311 V312 V313	5AQP 1V2 1V2 12AU7 6AQ5	CRT High Voltage Rectifier High Voltage Rectifier High Voltage Control Tube RF Oscillator	Adj.Vert.& Horiz.Gain(par.4-13,4-14 none none none none
V314	6DJ8/6BQ7A	a. +585-volt Series Regulator b. +585-volt Control Tube	none none

Sect. IV Page 8



4-8 ADJUSTMENT PROCEDURE

Usually a particular oscilloscope will not need complete testing and calibration. Only one or two tests will be needed and they can be done without completing the entire test procedure.

The following procedures are listed in a recommended sequence for a complete test and calibration operation. In general, tubes are the main cause of trouble and new ones should be tried before making adjustments or other component replacements.

Specifications for the \bigoplus Model 130B Oscilloscope are given in the front of this manual. The following test procedures contain extra checks to help you analyze a particular instrument. These extra checks and the data they contain can not be considered as specifications.

A fifteen minute warm-up and power supply output voltage measurements are always recommended before making any other test or adjustment.

4-9 TURN ON

When turning the oscilloscope on for the first time after repair in any circuit, measure resistance from power supplies to ground. They usually will be within 25% of the following:

110 ohms
0,000 ohms
9,000 ohms
5,000 ohms

CAUTION

When first turning an oscilloscope on after power supply repairs, turn the intensity and both positioning controls full counterclockwise before applying power. Failure to follow this precaution can cause permanent cathode-ray tube damage.

4-10 POWER SUPPLIES

The power supplies in the oscilloscope are extremely stable and will require infrequent adjustment. The output voltages should be measured at regular intervals but unnecessary adjustments should be avoided.

Power supply voltages may be measured at the points indicated in Figure 4-3.

To adjust the power supply section, refer to Figure 4-3, and proceed as follows:

A. LOW VOLTAGE SUPPLY

1) Turn sweep generator off by turning the HOR. SENS. switch to 50 volt/cm position.

2) Permit the 130B to warm up for at least five minutes at a line voltage of 115/230 volts.

3) Measure power supply voltages with line volts set to 115 volts. The voltages will normally be within the limits given in Table 4-3. Control R332 can be adjusted if necessary to set the +100 volt supply within limits.

If adjustment of the +100 volt supply was necessary, all sweep timing, calibrator and gain adjustments must be checked.

If poor low voltage supply regulation is suspected, the following check may be made:

- -- Check the regulation of each power supply voltage as the power line voltage is varied between 103 and 127 volts. All regulated voltages should remain within ±1% over this range of line voltage.
- Measure the ac ripple on each supply voltage. This ac voltage should not exceed the amount specified in Table 4-3.

Supply	Tolerance (115/230 volt line)	Variation ±109 line voltage change	% Nominal Ripple at 115/230V
+100 V	± 4%	±1%	5 mv
-150 V	± 4%	±1%	5 mv
+300 V	± 4%	±1%	5 mv
+585 V	± 4%	± 1%	60 mv

TABLE 4-3. REGULATED POWER SUPPLY TOLERANCES

If any output does not regulate or has excessive ripple, replace the Series Regulator Tube or the Control Tube of that supply. It must be kept in mind, however, that loss of regulation of the -150 volts will cause the other supplies to lose regulation, and that loss of regulation of the +100 volts will cause the +585 volt and +300 volt supplies to lose regulation also.
Sect. IV Page 10

SERVICING ETCHED CIRCUIT BOARDS

Excessive heat or pressure can lift the copper strip from the board. Avoid damage by using a low power soldering iron (50 watts maximum) and following these instructions. Copper that lifts off the board should be cemented in place with a quick drying acetate base cement having good electrical insulating properties.

Use only high quality rosin core solder when repairing etched circuit boards. NEVER USE PASTE FLUX. After soldering, clean off any excess flux and coat the repaired area with a high quality electrical varnish or lacquer.

A break in the copper should be repaired by soldering a short length of tinned copper wire across the break.

When replacing tube sockets it will be necessary to lift each pin slightly, working around the socket several times until it is free.



1. Apply heat sparingly to lead of part to be replaced. Remove part from card as iron heats the lead.



3. Bend clean tinned leads on new part and carefully insert through holes on board.



2. Using a small awl, carefully clean inside of hole left by old part.



4. Hold part against board and solder leads. Avoid overheating the board.

Figure 4-4

1) The -2550 volt output is measured on the resistor board under the base of the crt. Connect an appropriate dc voltmeter (such as an @ Model 410B VTVM with an @ Model 459A DC Resistive Voltage Multiplier) to the Junction of R338, C312 (marked -2550 CATH. on cover).

2) With the line at 115/230 volts the high voltage should measure $-2550 \pm 4\%$. Control R334 can be adjusted if necessary to set the -2550 volt supply within limits.

If poor -2550 voltage supply regulation is suspected the following check may be made:

- -- Check the regulation by varying the line voltage between 103 and 127 volts. The -2550 should remain within $\pm 1\%$ over this range of line voltage. If the -2550 supply does not regulate check the control tube V312.
- 3) Set SWEEP TIME/CM to 5 MILLISECONDS.
- 4) Set HOR. SENS. to INT. SWEEP X1.

5) Set SWEEP MODE fully clockwise to FREE RUN.

6) Set INTENSITY control to 10 o'clock.

7) Set Int. Adj., R343, until the trace is just visible.

8) Set SWEEP MODE fully counterclockwise to PRESET.

9) Set INTENSITY control for a low intensity spot.

10) Center spot and adjust FOCUS control and ASTIGMATISM (R303) to obtain a small round and sharply focused spot.

4-11 REPLACING AND ADJUSTING THE CRT

To replace the cathode-ray tube, refer to Figure 2-10, and proceed as follows:

1) Turn off and remove the 130B from the cabinet.

2) Loosen the clamp on the crt socket. (Cabinet model; remove cover from High Voltage terminal board to get access for screwdriver through board).

3) Remove the front-panel bezel.

4) With a screwdriver loosen the crt base from socket. Free the crt from the socket by pressing on the center of the tube base with one hand while supporting the front of the crt with the other.

DANGER - Do not apply force on neck of tube.

5) Remove the crt through the front panel.

<u>CAUTION</u> - HANDLE THE CATHODE RAY TUBE CAREFULLY.

6) Insert the replacement crt through the front panel and seat in socket.

7) Replace front-panel bezel.

8) Adjust the socket assembly so that the face of the crt just misses the bezel assembly. Tighten the clamp just enough to hold the crt in place loosely.

NOTE

Turn the INTENSITY control to minimum when first applying power to a crt. The phosphor can be damaged quickly by too much brightness.

9) Set the INTENSITY control fully counterclockwise. Turn the 130B on and allow to warm up.

10) Set the SWEEP MODE control to FREE RUN.

11) Adjust the INTENSITY control to obtain a weak trace; adjust the FOCUS control for a sharp trace, and with the vertical position control, center the trace vertically.

12) Align trace with graticule using the alignment handle at rear of crt.

<u>CAUTION</u> - Do not over-tighten crt clamp or tube damage may result.

13) Making certain the crt face is close to but not touching the bezel assembly, tighten the clamp on the crt socket only enough to hold the crt from turning. If the face of the tube touches the bezel assembly, Newton rings may be visible.

14) Readjust the astigmatism; see paragraph 4-10B.

15) Check the gain calibration of the Vertical and Horizontal Amplifiers by setting the VERTICAL and HORIZONTAL SENSITIVITY selectors to CAL, and if necessary, adjusting R40 (Figure 4-5) to obtain 6 cm vertical deflection and R144 (Figure 4-6) to obtain 6 cm horizontal deflection on the trace; see paragraph 4-13B and 4-14B. Sect. IV Page 12



1) Set HORIZ. SENSITIVITY to INT. SWEEP X1; SWEEP TIME/CM to 1 MILLISECOND; SWEEP MODE to PRESET; SYNC. to INT.; and TRIG-GER LEVEL to "0".

NOTE

If PRESET (SWEEP MODE) is maladjusted, you may not obtain a trace. See paragraph 4-16.

2) Set the VERT. SENSITIVITY switch to 50 MILLIVOLTS/CM. Place the VERNIER control in CAL.

3) Connect the signal source to the vertical IN-PUT terminals and set its rms output voltage, read on the vtvm, to 106 millivolts (300 mv/ $2\sqrt{2}$) and its output frequency to 1000 cps.

4) Adjust R40 (see Figure 4-5) to obtain exactly 6 centimeters deflection.

5) Set the VERT. SENSITIVITY switch to CAL.

6) Adjust the R240 (see Figure 4-7) for exactly 6 centimeters deflection.

4-13 ADJUSTING THE VERTICAL AMPLIFIER

The following adjustments are located in the vicinity of the Vertical Amplifier or the VERT. SENSITIVITY switch as shown in Figure 4-5.

A. VERNIER BALANCE ADJUSTMENT

To adjust VERNIER balance, allow the instrument to warm up 15 minutes and adjust Vertical balance as shown by Figure 2-2. Then refer to Figure 4-5 and proceed as follows:

1) Short the INPUT terminals and set the INPUT switch to DC.

2) Set VERT. SENSITIVITY to 1 MILLIVOLT/CM, and VERNIER to CAL.

3) Center spot (or trace) with VERT, POS, control.

4) Turn VERNIER fully counterclockwise and return spot to center with R20, the Bal. Adj.

The trace should now be stationary as the VER-NIER is rotated.

B. VERT. AMPL. GAIN AND FREQ. RESP. AD-JUSTMENTS

To adjust the Vertical Amplifier gain and frequency response adjustments refer to Figure 4-5 and proceed as follows:

1) Adjust Vertical and VERNIER balance as indicated by Figure 2-2 and paragraph 4-13A.

2) Set HORIZ. SENSITIVITY to INT. SWEEP X1 and set the SWEEP TIME/CM switch to 2 MILLI-SECONDS, SWEEP MODE to PRESET and TRIG-GER LEVEL to "0".

3) Place VERTICAL SENSITIVITY switch and its VERNIER in CAL.

4) Adjust R40 for exactly 6 centimeters deflection.

5) Set VERT. SENSITIVITY to 50 MILLIVOLTS/ CM. Set SYNC to INT.

6) Set SWEEP TIME/CM switch to 5 MICRO-SECONDS.

7) Connect a 50 kc square wave to the Vertical INPUT and adjust the square-wave amplitude for 6 to 8 centimeters deflection.

8) Adjust C12 for best square wave. To give C12 a range of adjustment sufficient to compensate for variations of tube characteristics, C13 may be connected in parallel with C12 to increase the maximum capacity to $1340 \mu\mu$ f.

C. INPUT ATTENUATOR FREQUENCY RESPONSE ADJUSTMENTS

To adjust frequency response of the input attenuator refer to Figure 4-5 and proceed as follows:

1) Connect a 5 kc square wave to the Vertical INPUT.

2) Set SWEEP TIME/CM to obtain 3 or 4 cycles of the square wave.

3) Make the indicated adjustment for best squarewave presentation on the following ranges:

VOLT/CM	ADJUST
10	C2
1	C4
.1	C3



4-14 ADJUSTING HORIZONTAL AMPLIFIER

To adjust the Horizontal Amplifier, refer to Figure 4-6 and proceed as follows:

A. VERNIER BALANCE ADJUSTMENT

To adjust the VERNIER balance allow the instrument to warm up thoroughly and adjust Horizontal balance as shown by Figure 2-3; then refer to Figure 4-6 and:

1) Short the INPUT terminals and set the INPUT switch to DC.

2) Set the HORIZ. SENSITIVITY to 1 MILLI-VOLT/CM and the VERNIER to CAL.

3) Center the spot with the HORIZ. POS. control.

4) Turn the VERNIER completely counterclockwise and return the spot to the center with R120, Bal. Adj. The spot will now be stationary as the VERNIER is rotated.

B. AMPLIFIER GAIN AND FREQUENCY RESPONSE ADJUSTMENTS

To adjust the gain and frequency response, refer to Figure 4-6 and proceed as follows:

1) Set HORIZ. SENSITIVITY switch to CAL. and the VERNIER to CAL.

2) Adjust R144 for exactly 6 centimeters deflection.

3) Set VERT. SENSITIVITY to 2 VOLTS/CM.

4) Connect an 8 kc (approximately) sine wave to the Vertical INPUT of the oscilloscope and to the SYNC. IN terminal of the 211A square wave generator; adjust the sine wave for 10 cm deflection.

5) Set HORIZ. SENSITIVITY to 50 MILLI-VOLTS/CM.

6) Connect a 50 kc square wave to the Horizontal INPUT, and adjust the square wave amplitude for 6 to 8 cm deflection.

7) Adjust C114 for best square wave response.

NOTE

Some vacuum tubes require more capacity for compensation than the maximum value of Cl14.

Capacitor C115 may be connected in parallel with C114 to increase the maximum capacity to 1340 $\mu\mu$ f, permitting a greater percentage of vacuum tubes to be used.

C. INPUT ATTENUATOR FREQUENCY RESPONSE ADJUSTMENTS

To adjust the frequency response at the input attenuator, refer to Figure 4-6 and proceed as follows:

1) Set VERT. SENSITIVITY to 2 VOLTS/CM.

2) Connect an 800 cps (approximately) sine wave to the Vertical INPUT of the oscilloscope and to the Sync-In terminal of the 211A square wave generator; adjust the sine wave for 10 centimeter deflection.

3) Set HORIZ. SENSITIVITY to 10 VOLTS/CM (VERNIER in CAL.).

4) Connect a 5 kc square wave to the Horizontal INPUT and adjust its amplitude for 6 centimeter deflection.

5) Make the following adjustments on the ranges indicated for the best square wave response, adjusting the square-wave amplitude to 6 centimeters on each range.

ADJUST (see Fig. 4-6)
C102
C104
C103

4-15 PHASE SHIFT ADJUST

Phase shift between Vertical and Horizontal Amplifiers.

If the square wave response of the Vertical and Horizontal Amplifier was carefully set, the relative phase shift between the two amplifiers should not exceed one degree at frequencies below 50 kc.

To check Phase Balance:

1) Set VERT. and HORIZ. SENSITIVITY to 50 MILLIVOLTS/CM and VERNIER to CAL.

2) Apply a 50 kc sine-wave signal to the HORIZ. INPUT and VERT. INPUT. Center pattern and adjust signal amplitude for 6 cm vertical and 6 cm horizontal deflection.





Figure 4-7. Sweep Generator Adjustment Locations

The opening of the pattern, if any, should not exceed a tenth of a centimeter. If necessary to correct phase difference, adjust C114 for closure of the pattern (Figure 4-6).

4-16 ADJUSTING PRESET

To adjust Preset refer to Figure 4-7 and proceed as follows:

1) Set VERT. SENSITIVITY to OFF, SWEEP TIME/CM switch to .1 MILLISECOND and SWEEP MODE control maximum counterclockwise to PRE-SET. Set SYNC selector to INT.

2) Connect a dc voltmeter between ground (+) and the center tap of R220 the PRESET adjust control.

3) Turn R220 fully counterclockwise. Then slowly adjust R220 clockwise until the sweep generator begins to free fun. Turn R220 counterclockwise until the sweep just stops and record this voltage which should be about -26 volts.

4) Set R220 to give a voltmeter indication exactly 2 volts less negative than the voltage noted.

4-17 ADJUSTING THE SAWTOOTH GENERATOR AND SWEEP AMPLIFIER

To adjust the Sawtooth Generator and Sweep Amplifier refer to Figures 4-6 and 4-7 and proceed as follows:

1) Set SYNC to INT., SWEEP MODE to PRESET and TRIGGER LEVEL to "0".

2) Set HORIZ. SENSITIVITY to INT. SWEEP X1, SWEEP TIME/CM to 1 MILLISECOND, and its VERNIER to CAL.

3) Connect 1 kc (1000 μ sec) time markers to the Vertical INPUT.

4) Set R261, 1 Millisecond Adj., on the potentiometer board, to its mechanical center.

5) Adjust R134, Sweep Attenuator, (Figure 4-6) and HORIZ. POS. for approximately 1 time marker/ cm. This is a rather coarse adjustment. Set it as close as is practical. Then make the final adjustment with R261.

6) Set SWEEP TIME/CM to 10 MILLISECONDS and connect 100 cycle (10,000 μ sec) time markers to the Vertical INPUT.

7) Adjust R260 (Figure 4-7) for 1 marker per centimeter.

8) Set SWEEP TIME/CM to .1 SECOND and connect 10 cycle (100,000 μ sec) time markers to the Vertical INPUT.

9) Adjust R259 (Figure 4-7) for 1 marker per centimeter.

10) Disconnect the time mark generator from the Vertical INPUT, set VERT. SENSITIVITY switch to 10 VOLTS/CM, the Vertical and Horizontal input switches to AC, and SWEEP TIME/ CM to 5 MICROSECONDS.

11) Disconnect the blue-white lead (Figure 4-7) from the Sweep Generator board and connect it through a 1 microfarad capacitor to the Horizontal INPUT.

12) Connect a wire between V206 pin 1 (6AW8) and the Vertical INPUT.

13) Connect a 50 kc square wave from the 600 ohm output of the @ Model 211A to the Horizontal INPUT and adjust its amplitude for about 6 centimeter deflection.

14) Adjust C113 near the HORIZ. SENSITIVITY switch, and C116 and C118 on the Horizontal Amplifier board (Figures 4-6 and 4-7), for best square wave presentation. Remove the wire between V206 pin 1 and the Vertical INPUT.

15) Reconnect the blue-white lead to the Sweep Generator board, and connect the time marker generator to the Vertical INPUT. Make the following adjustment as indicated for 1 time marker per centimeter.

Time Marker	SWEEP TIME/CM	(Fig. 4-7) Adjust
1 µsec (1 mc)	1 MICROSEC.	C227
10 µsec (100 kc)	10 MICROSEC.	C225
.1 msec(10 kc)	.1 MILLISEC.	C223

16) Set HORIZ. SENSITIVITY to INT. SWEEP X5, and adjust R164 (Figure 4-6), for markers 5 centimeters apart.

17) Connect a 500 kc signal to the Horizontal INPUT, set SWEEP TIME/CM to 1 MILLISECOND, set SYNC to EXT, and adjust R229, Sweep Length, (Figure 4-7) for a trace about 10.5 centimeters long.

SCHEMATIC DIAGRAM NOTES

- 1. Heavy solid line shows main signal path; heavy dashed line shows control, secondary signal, or feedback path.
- 2. Heavy box indicates front-panel engraving; light box indicates chassis marking.
- 3. Arrows on potentiometers indicate clockwise rotation as viewed from the round shaft end, counterclockwise from the rectangular shaft end.
- 4. Resistance values in ohms, inductance in microhenries, and capacitance in micromicrofarads unless otherwise specified.
- 5. Rotary switch schematics are electrical representations; for exact switching details refer to the switch assembly drawings.
- 6. Relays shown in condition prevailing during normal instrument operation.
- 7. ‡ indicates a selected part. See parts list.
- 8. Interconnecting parts and assemblies are shown on cable diagram.
- 9. * indicated value adjusted at factory. Part may be omitted.

VOLTAGE AND RESISTANCE DIAGRAM NOTES

1. Each tube socket terminal is numbered and lettered to indicate the tube element and pin number, as follows:

*	=	no tube element	Р	=	plate
H	=	heater	Т	=	target (plate)
K	=	cathode	R	=	reflector or repeller
G	=	control grid	Α	=	anode (plate)
Sc	=	screen grid	S	=	spade
Sp	=	suppressor grid	Sh	=	shield
Hm	=	heater mid-tap	NC	=	no external connection to socket
IS	=	internal shield	Δ	=	indefinite reading due to circuit (see 2.)

The numerical subscript to tube-element designators indicates the section of a multiple-section tube; the letter subscript to tube-element designators indicates the functional difference between like elements in the same tube section, such as t for triode and p for pentode.

A socket terminal with an asterisk may be used as a tie point and may have a voltage and resistance shown.

- 2. Voltages values shown are for guidance; values may vary from those shown due to tube aging or normal differences between instruments. Resistance values may vary considerably from those shown when the circuit contains potentiometers, crystal diodes, or electrolytic capacitors.
- 3. Voltage measured at the terminal is shown above the line, resistance below the line; measurements made with an electronic multimeter, from terminal to chassis ground unless otherwise noted.
- 4. A solid line between socket terminals indicates a connection external to the tube between the terminals; a dotted line between terminals indicates a connection inside the tube. Voltage and resistance are given at only one of the two joined terminals.







S2 Q VERT. SENSITIVITY

HORIZONTAL AMPLIFIER VOLTAGE - RESISTANCE DIAGRAM (VIEWED FROM ETCHED SIDE)

V101 (12AU7) V104 (6807/6DJ8) DIFFERENTIAL AMPLIFIER INVERTER AMPLIFIER 420 140K +1.8 V ov 00 120 K +160V, 0VAC +248 V +54V 50K NC 140K +160V, 6.3VAC +252 + 64 V 65 16 K +252 V NC + 51V +54 V +248V 140 K 0 +1.8V 120K 100K +420V 140K V104 V103 V102 V101) LD-E-118 LD-S - 189 1 1 V102 (12AU7) V103 (12AT7) DIFFERENTIAL AMPLIFIER DIFFERENTIAL AMPLIFIER +138V 90K + 53V 140K +56V 7.5K +24PV +140V 50K +158V NC 35K +88V 90 +82V 90 +76V 80 H +76V + 158V +64V 65 +140V +248V K2 G2 G2 Ka 50K +53V +56V +138V 2200 LD-E-1308CD FIGURE 4-10



CI12, CI17, C124 RI13, R124, R126, R127, R135, R156-163 SWEEP GENERATOR VOLTAGE - RESISTANCE DIAGRAM (VIEWED FROM RIGHT SIDE)





V202 (12AT7)



V203 (6U8)

SWEEP START-STOP TRIGGER



V204 (6C4)

GATE OUT CATHODE FOLLOWER







V205 (12AL5)





V201 (6BQ7A/6DJ8) TRIGGER AMPLIFIER

V207 (12AX7) RETRIGGERING BIAS CONTROL







V206 (6AW8)

INTEGRATOR CATHODE FOLLOWER



C112, C117, C124 R113, R124, R126, R127, R135, R156-163





This drawing is intended for the operation and maintenance of Hewlett-Packard equipment and is not to be used otherwise cr reproduced withan consent of the

reproduced without written consent of the Hewlett-Packard Company. 1308 - ST/CM SWITCH - T946A SWEEP GENERATOR VOLTAGE - RESISTANCE DIAGRAM (VIEWED FROM RIGHT SIDE)



V202 (12AT7) TRIGGER GENERATOR V203 (608) SWEEP START-STOP TRIGGER V204 (6C4) GATE OUT CATHODE FOLLOWER





1 1

V203 V204 V205

V207

V206

I

1

1

V202

V201



V205 (12AL5) INTEGRATOR SWITCH

> K₂ ⊢NC

+88V

+98V

6.1

90



| | | V201 (6BQ7A/6DJ8)







V207 (12AX7)

RETRIGGERING BIAS CONTROL



V206 (6AWB) INTEGRATOR CATHODE FOLLOWER

+ 100 V

98V

120K

+ 98V

800

FIGURE 4-12









CATHODE RAY TUBE WARRANTY

The cathode ray tube supplied in your Hewlett-Packard Oscilloscope and replacement cathode ray tubes purchased from \bigoplus , are guaranteed against electrical failure for one year from the date of sale by the Hewlett-Packard Company. Broken tubes or tubes with burned phosphor are not included in this guarantee.

Your local Hewlett-Packard representative maintains a stock of replacement tubes and will be glad to process your warranty claim for you. Please consult him.

Whenever a tube is returned for a warranty claim, the reverse side of this sheet must be filled out in full and returned with the tube. Follow shipping instructions carefully to insure safe arrival, since no credit can be allowed on broken tubes.

SHIPPING INSTRUCTIONS

- Carefully wrap the tube in 1/4" thick cotton batting or other soft padding material.
- 2) Wrap the above in heavy kraft paper.
- 3) Pack in a rigid container which is at least 4 inches larger than the tube in each dimension.
- 4) Surround the tube with at least four inches of packed excelsior or similar shock absorbing material. Be certain that the packing is tight all around the tube.
- 5) Tubes returned from outside the continental United States should be packed in a wooden box.
- Ship prepaid preferably by AIR FREIGHT or RAILWAY EXPRESS. We do not recommend parcel post or air parcel post shipment.

HEWLETT-PACKARD CO. 275 PAGE MILL ROAD, PALO ALTO, CALIF. U.S.A.

	CRT WARRANTY CLAIM
FROM:	DATE:
NAME:	A ABAR / CONTRACTOR
COMPANY:	and the second sec
ADDRESS:	
Person to contact for further info	mation:
TITLE:	
COMPANY:	
ADDRESS:	
	1-113-007-74
To process your claim quickly ple	ase enter the information indicated below:
1) @ INSTRUMENT MODEL	SERIAL
2) TUBE TYPE	SERIAL
3) ORIGINAL TUBE	REPLACEMENT TUBE
4) YOUR PURCHASE ORDER NO	
5) DATE PURCHASED	
6) PURCHASED FROM	
7) COMPLAINT: (Please describ	e nature of trouble)
	and the set of the set of the second
B) OPERATING CONDITIONS: (F	lease describe conditions prior to and at time of failure
	SIGNATURE
	SIGNATURE

SECTION V TABLE OF REPLACEABLE PARTS

Readily available standard-components have been used in this instrument, whenever possible. However, special components may be obtained from your local Hewlett-Packard representative or from the factory.

NOTE

When ordering parts always include:

- 1. 👜 Stock Number.
- 2. Complete description of part including circuit reference.
- Model number and serial number of instrument.
- 4. If part is not listed give complete description, function, and location of part.

If there are any corrections for the Table of Replaceable Parts they will be listed on an Instruction Manual Change sheet at the front of this manual.

Sect.V Page 2

Model 130B

TABLE OF REPLACEABLE PARTS

REF.	DESCRIPTION, MFR. * & MFR. DESIGNATION	STOCK NO.	#		
C1	Capacitor: fixed, mylar, 0.1 μ f ±20%, 600 vdcw Texas Capacitor Co.	16-110	2		
C2	Capacitor: variable, ceramic, 5-20 μμf 500 vdcw L*	13-20	2		
C3	Capacitor: variable, mica, 14-50 μμf 500 vdcw L*	13-33	4		
C4	Capacitor: variable, mica, 50-380 μμf 175 vdcw QQ*	13-5	2		
C5	Capacitor: fixed, silver mica, 820 μμf ±5%, 500 vdcw Z*	15-104	4	ŝ	
C6	Capacitor: fixed, mica .01 μ f ±5%, 300 vdcw Z*	14-24	3	1	
C7, 8	Capacitor: fixed, ceramic, 10,000 $\mu\mu f \pm 20\%$, 1000 vdcw CC*	15-43	6		
C9	Capacitor: fixed, mica, 75 μμf ±5%, 300 vdcw V*	14-75	3		
C10	Capacitor: fixed, mica, 15 $\mu\mu f \pm 10\%$, 500 vdcw V*	14-15	2		
C11	This circuit reference not assigned				
C12	Capacitor: variable, mica 170-780 µµf 175 vdcw QQ*	13-32	3		
C13	Capacitor: fixed, mica, 560 μμf ±10%, 500 vdcw Electrical value adjusted at factory	14-81	2		
C14	This circuit reference not assigned		. 1	10	
C15, 16	Capacitor: fixed, titanium dioxide, 2 $\mu\mu$ f ±5%, 500 vdcw Electrical value adjusted at factory DD*	15-118	4		
C17, 18	Capacitor: fixed, titanium dioxide, 3.3 $\mu\mu$ f ±10%, 500 vdcw DD*	15-78	4		
C19	Capacitor: fixed, titanium dioxide, 1.5 $\mu\mu$ f ±20%, 500 vdcw DD*	15-38	4		
C20	This circuit reference not assigned				

CIRCUIT REF.	DESCRIPTION, MFR. * & MFR. DESIGNATION	STOCK NO.	#	(Pic		
C21	Same as C19 Electrical value adjusted at factory					
C22	Same as C7	1.634350				
C23 thru C100	These circuit references not assigned					1. j
C101	Same as C1				1.1	W. al
C102	Same as C2	1.000		1.7		
C103	Same as C3			1	1.1	1.1.1
C104	Same as C4					1.00
C105	Same as C5	14-22		1.11	F	
C106	Same as C6	New York		17.14	. 321	
C107, 108	Same as C7					
C109	Same as C9					
C110	Same as C10					5 240
C111	Capacitor: fixed, mica, 200 μ f ±10%, 500 vdcw V*	14-200	1			
C112	This circuit reference not assigned	1.244		100		1.2
C113	Capacitor: variable, ceramic, 7-45 μμf, 500 vdcw L*	13-1	3			
C114	Same as C12					1
C115	Same as C13					
C116	Same as C3	1.5253			15.1	1.500
C117	This circuit reference not assigned			1	inst-	
C118	Same as C3	1		14A		i
C119, 120	Same as C19 Electrical value adjusted at factory					33
C121, 122	Same as C17				1 12	N Agent
Land La	10-1 30 38 L 34	and and and and		i		22

TABLE OF REPLACEABLE PARTS

TABLE OF REPLACEABLE PARTS

DE	D	DESCRIPTION, MFR. * & MFR. DESIGNATIO	ON	STOCK NO.	#				
apa 2. El	Ca	pacitor: fixed, titanium dioxide, 2.2 $\mu\mu$ f ±10%, 500 vdcw Electrical value adjusted at factory D	D*	15-52	2				
nis	Th	is circuit reference not assigned							
me El	Sai	ne as C123 Electrical value adjusted at factory		in an		112			
apa . 0	Ca	pacitor: fixed, paper, .01 μf ±10%, 600 vdcw C	C*	16-11	2				
nes	Th	ese circuit references not assigned				an I	'n,		
apa . C	Ca	pacitor: fixed, paper, .047 μf ±10%, 1000 vdcw C	C*	16-114	2				
nis	Th	is circuit reference not assigned		8.1 U			100	1	
me	Sai	ne as C201				-	1	3	
apa 27	Ca	pacitor: fixed, mica, 27 $\mu\mu f \pm 10\%$, 500 vdcw	v*	14-17	2				
apa 68	Ca	pacitor: fixed, mica, 68 $\mu\mu f \pm 10\%$, 500 vdcw	v*	14-60	1				
me	San	me as C7		22.18					
apa 10	Ca	pacitor: fixed, ceramic, 10 μμf ±0.5%, 500 vdcw	K*	15-30	1		1		
me	Sai	me as C204							
apa 20	Ca	pacitor: fixed, silver mica, 200 µf ±5%, 500 vdcw	Z*	15-103	1				
apa . 0	Ca	pacitor: fixed, paper, .0022µf ±10%, 600 vdcw C	C*	16-22	2				
apa . 0	Ca	pacitor: fixed, ceramic, .02 μf -0% +100%, 600 vdcw	G*	15-85	2				
me	Sai	me as C7		· · · · · · · · · · · · · · · · · · ·					
me	Sai	me as C210				-			
me	Sai	ne as C15		2					
ime	Sai	ne as C15							

ERRATA:

C218:

.8: Add capacitor, fixed, ceramic, .Oluf +20%, e5 1000 vdcw; -hp- Stock No. 15-43, Mfr., CC

IABLE OF REPLACEABLE PARTS

CIRCUIT REF.	DESCRIPTION, MFR. * & MFR. DESIGNATION	STOCK NO.	#		
C216	Capacitor: fixed, ceramic, .01 μf ±20%, 250 vdcw CW*	15-135	1		
C217 thru C219	These circuit references not assigned				
C220	Capacitor: fixed, mylar, 1.0 μ f ±5%, 200 vdcw CW*	16-102	2		
C221	Capacitor: fixed, mylar, 0.1 μf ±5%, 200 vdcw CW*	16-103	1		
C222	Capacitor: fixed, mylar, .01 μf ±5%, 400 vdcw CW*	16-101	1		
C223	Same as C4				
C224	Same as C5				
C225	Same as C113				
C226	Same as C9				
C227	Same as C113				
C228	Same as C220				
C229	Capacitor: fixed, paper, 0.1 μf ±10%, 400 vdcw CC*	16 -35	4		
C230	Same as C126				
C231	Capacitor: fixed, mica, 680 μμf ±10%, 500 vdcw Z*	14-21	1		
C232	Capacitor: fixed, mica, 470 $\mu\mu$ f ±10%, 500 vdcw V*	14-62	1		
C233 thru C300	These circuit references not assigned				
C301	Capacitor: fixed, electrolytic, 2 section, 20 μ f/sect.CC*	18- 22 HP	2		
C302	Capacitor: fixed, electrolytic, 80 μf, 300 vdcw X*	18-64	1		
C303	Capacitor: fixed, electrolytic, 160 μ f, 450 vdcw Electrical value adjusted at factory CC*	18-51HP	1		

* See "List of Manufacturers Code Letters For Replaceable Parts Table".

Total quantity used in the instrument.

CIRCUIT REF.	DESCRIPTION, MFR. * & MFR. DESIGNAT	ION	STOCK NO.	#			
C304A, B	Same as C301						
C305	Capacitor: fixed, paper, 0.1 μ f ±10%, 1000 vdcw	J*	16-82	1			
C306, 307	Capacitor: fixed, tubular, .01 μ f ±20%, 400 vdcw	cc*	16-109	5			
C308	Same as C229						
C309	Same as C306						
C310A, B	Capacitor: fixed, electrolytic, 4 sections, 20 μ f/sect., 450 vdcw	cc*	18-42HP	2			
C311, 312	Capacitor: fixed, paper, .0068 μ f ±10%, 5000 vdcw	cc*	16-93	3			
C313	Capacitor: fixed, paper, .015 μ f ±10%, 3000 vdcw	cc*	16-126	1			
C314, 315	Capacitor: fixed, paper, .0015 μ f ±20%, 500 vdcw	cc*	16-125	2			
C316	Capacitor: fixed, ceramic, .02 μ f -0% + 100%, 600 vdcw Radio Mat.	Co.	15-85	1	100		
C317 thru							
C319	These circuit references not assigned						
C320	Same as C306						
C321	This circuit reference not assigned						
C322	Capacitor: fixed, ceramic, 2000 $\mu\mu$ f \pm 20%, 1000 vdcw	L*	15-80	2			
C323	Capacitor: fixed, mica, 200 $\mu\mu$ f ±10%, 500 vdcw	v*	14-200	1			
C324	Same as C229						
C325	Same as C306						
C326	Capacitor: fixed, paper, .001 μ f ±10%, 600 vdcw	cc*	16-21	1			
C327	Same as C7	¢					
C 32 8	Capacitor: fixed, paper, .22 μ f 400 vdcw	cc*	16-48	1			
							, I

TABLE OF REPLACEABLE PARTS

CIRCUIT REF.	DESCRIPTION, MFR. * & MFR. DESIGNATIO	ON	STOCK NO.	#		
CR301, 302	Rectifier, selenium	N*	212-114	2		
CR303 thru CR306	Rectifier, silicon: 500 V PIV B	3V*	212-147	4		
F301	Fuse, cartridge: 2 amp, slow blow for 115 V 1/2 amp, slow blow, for 230 V operation	Т*	211-16 211-20	1 1		
11, 2	Lamp, glow: 1/25 W, 90 vdcw, 65 VAC, NE2	0*	211-43	6		
I3 thru I100	These circuit references not assigned					
1101, 102	Same as I1					
I103 thru I301	These circuit references not assigned					
I302 thru I305	Lamp, incandescent: 6-8V, .15 amp	0*	211-47	4		
1306	This circuit reference not assigned					
1307, 308	Same as I1					
J1	Binding Post Assembly: red H	[P *	AC-10D	4		
	Insulator, binding post (rack model) H	P *	AC-54A	4		
	Insulator, binding post (cabinet model) H	(P *	AC-54B	2		
	Connector Assembly H	[P *	G-76J	1		
J2	Connector, receptacle: male, 3 contact (on rear panel, rack mount only) C	'U*	125-58	2		
J3 thru J100	These circuit references not assigned					

TABLE OF REPLACEABLE PARTS

Sect. V Page 8

CIRCUIT REF.	DESCRIPTION, MFR. * & MFR. DESIGNA	TION	STOCK NO.	#			
J101	Binding Post Assembly: red	HP*	AC-10D	(2)	See	J1	
	Insulator, binding post (rack model)	HP*	AC-54A	(2)	See	J1	
	Insulator, binding post (cabinet model)	HP*	AC-54B	(1)	See	J1	
	Connector Assembly	HP*	G-76K	1			
J102	Same as J2						
L201	Coil, choke: $1 \ \mu h \pm 10\%$	CG*	48-96	1			
L202	Coil, r.f.: 360 µh	CG*	48-63	1			
L203 thru L300	These circuit references not assigned						
L301	Coil, r.f.: 5 µh	CG*	48-61	1			
L302	Coil, r.f.: 200 μ h	CG*	48-34	1			
P1	Cord, power Cornish V	Wire Co	812-106	1			
R1	Resistor: fixed, deposited carbon, 900,000 ohms ±1%, 1 W	NN*	33-900K	2			
R2	Resistor: fixed, deposited carbon, 90,000 ohms $\pm 1\%$, 1/2 W	NN*	33-90K	2			
R3	Resistor: fixed, deposited carbon, 9000 ohms $\pm 1\%$, 1/2 W	NN*	33-9000	1			
R4	Resistor: fixed, deposited carbon, 1000 ohms $\pm 1\%$, 1/2 W	NN*	33-1000	2			
R5	Resistor: fixed, deposited carbon, 1 megohm $\pm 1\%$, 1/2 W	NN*	33-1M	4			
R6, 7	Resistor: fixed, composition, 100,000 ohms $\pm 10\%$, 1/2 W	в*	2 3-100K	4			
R8, 9	Resistor: fixed, composition, 100 ohms $\pm 10\%$, 1/2 W	В*	23-100	17			

TABLE OF REPLACEABLE PARTS

* See "List of Manufacturers Code Letters For Replaceable Parts Table". # Total quantity used in the instrument.

Model 130B

CIRCUIT REF.	DESCRIPTION, MFR. * & MFR. DESIGNATION	Designment of the stock no.	#			
R10A, B	Resistor: variable, dual concentric, linear taper, Rear section: 250 ohms ±10% Front section: 10,000 ohms ±10%, 2 W BO*	210-199	2			
R11	Resistor: fixed, composition, 2.7 megohms ±10%, 1/2 W B*	23-2.7M	3	1		
R12	Resistor: fixed, deposited carbon, 75,000 ohms ±1%, 1 W NN*	31-75K	2		off	
R13	This circuit reference not assigned			6		
R14, 15	Resistor: fixed, metal film, 100,000 ohms ±1%, 1/2 W	330-100K	4			
	The Daven Co.					
R16, 17	Resistor: fixed, deposited carbon, 30,500 ohms ±1%, 1 W NN*	31-30. 5K	4		1	
R18	Resistor: variable, composition, 20,000 ohms Part of S3 BO*	210-195	1	2.2		
R19	Resistor: fixed, deposited carbon, 15,000 ohms ±1%, 1/2 W NN*	33-15K	8			
R20	Resistor: variable, linear taper, 500 ohms $\pm 20\%$, 2/10 WBO*	210-203	2			
R21	Same as R19	ne dan m		1485-01		
R22	Resistor: fixed, deposited carbon, 1200 ohms ±1%, 1/2 W NN*	33-1200	2			
R23	Resistor: fixed, deposited carbon, 1800 ohms ±1%, 1/2 W NN*	33-1800	2			
R24	This circuit reference not assigned					
R25	Resistor: fixed, deposited carbon, 3000 ohms $\pm 1\%$, $1/2$ W NN*	33-3000	3			
R26, 27	These circuit references not assigned					
R28	Resistor: fixed, deposited carbon, 6000 ohms $\pm 1\%$, $1/2$ WNN*	33-6000	2			
R29, 3 0	Same as R3					
R31, 32	Same as R19					
		and the second states of			-	

TABLE OF REPLACEABLE PARTS

R33, 34 R35, 36 R37	Same as R8 Resistor: fixed, deposited carbon, 41,500 ohms $\pm 1\%$, 1 W Resistor: fixed, deposited carbon, 49,000 ohms $\pm 1\%$, 1 W	NN*	31-41. 5K			10	
R35, 36 R37	 Resistor: fixed, deposited carbon, 41,500 ohms ±1%, 1 W Resistor: fixed, deposited carbon, 49,000 ohms ±1%, 1 W 	NN*	31-41.5K				1
R37	Resistor: fixed, deposited carbon, 49,000 ohms $\pm 1\%$, 1 W			4			
		NN*	31-49К	2			
R38	Resistor: variable, composition, linear taper, 5000 ohms	G*	210-15	2			
R39	Resistor: fixed, composition, 560 ohms ±10%, 1/2 W	В*	23-560	2			
R40	Resistor: variable, composition, 5000 ohms ±30%, 1/3 W, linear taper	BO*	210-134	2			
R41	Resistor: fixed, deposited carbon, 800,000 ohms $\pm 1\%$, 1/2 W	NN*	33-800K	2			
R42, 43	Same as R8		a dan sadata		1.12	ň., 1	
R44	Same as R41				1		
R45	Resistor: fixed, composition, 5100 ohms ±5%, 1 W	в*	24-5100-5	1			
R46	Resistor: fixed, composition, 22,000 ohms $\pm 10\%$, 2 W	в*	25-22К	2			
R47	Resistor: fixed, composition, 10,000 ohms ±10%, 2 W	в*	25-10K	6	2.0		
R48	Same as R46					SIK.	14
849	Same as R45						
850, 51	Same as R8				- 7		
R52 thru R100	These circuit references not assigned				1		144
R101	Same as R1						a an
R102	Same as R2						
R103	Resistor: fixed, deposited carbon, 9000 ohms $\pm 1\%$, 1/2 W	NN*	33-9000	3	8		e i se

TABLE OF REPLACEABLE PARTS

CIRCUIT REF.	DESCRIPTION, MFR. * & MFR. DESIGNATION	STOCK NO.	#		0	
R104	Same as R4	an in a sid		-	84 J.	
R105	Same as R5	bela sole			101	
R106, 107	Same as R6	1			100	
R108, 109	Same as R8		1		1 Inte	
R110A,B	Same as R10					
R111	Same as R11					
R112	Same as R12	· · · · · · · · · · · · · · · · · · ·				1
R113	This circuit reference not assigned				2	1
R114, 115	Same as R14	and and a second	1		Lab 1	1.010
R116, 117	Same as R16	N - 10 - 1			1	1.1
R118	Resistor: variable, composition, 20,000 ohms BO*	210-195	1			
R119	Same as R19	ana alifait a	A DRU	-		
R120	Same as R20		dete			
R121	Same as R19	and the second		123	ene pro	elena:
R122	Same as R22	a financia de la composición de la comp				1
R123	Same as R23	a - 1911	failes.	AHE. 1		
R124	This circuit reference not assigned			010		0.00
R125	Same as R25					20
R126, 127	These circuit references not assigned				*	
R128	Same as R28			0.54		
R129, 130	Same as R103			1		
R131, 132	Same as R19					11.3
R133	Resistor: fixed, deposited carbon, 50,000 ohms ±1%, 1/2 W NN*	33-50K	1			and the second
R134	Resistor: variable, linear taper, 125,000 ohms BO*	210-110	1			
					1.11	

TABLE OF REPLACEABLE PARTS

Sect. V Page 12

CIRCUIT REF.	DESCRIPTION, MFR. * & MFR. DESIGNA	ATION	STOCK NO.	#		-	_
R135	This circuit reference not assigned						
R136	Resistor: fixed, deposited carbon, 900,000 ohms $\pm 1\%$, 1/2 W	NN*	33-900K	1			
R137, 138	Same as R8						
R139, 140	Same as R35						
R141	Same as R37						
R142	Same as R38						
R143	Same as R39						
R144	Same as R40						
R145	Resistor: fixed, deposited carbon, 1.5 megohms $\pm 1\%$, 1/2 W	NN*	33-1.5M	3			
R146, 147	Same as R8	-here	anoy oli				
R148	Same as R145						
R149	Resistor: fixed, composition, 6200 ohms ±5%, 2 W	в*	25-6.2K-5	1	[]===_]		
R150, 151	Resistor: fixed, composition, 7500 ohms $\pm 5\%$, 1 W	в*	24-7.5K-5	2			
R152, 153	Resistor: fixed, composition, 27,000 ohms ±10%, 2 W	в*	25-27K	2			
R154, 155	Same as R8					} (
R156 thru R163	These circuit references not assigned						
R164	Resistor: variable, composition, 2500 ohms ±20%, 1/2 W	G*	2 10-98	1			
R165	Resistor: fixed, composition, 56 megohms $\pm 10\%$, 1/2 W	в*	23-56M	1			
R166 thru R200	These circuit references not assigned		1.2				
			- 10 ²				

TABLE OF REPLACEABLE PARTS

Total quantity used in the instrument.

Т

CIRCUIT REF.	DESCRIPTION, MFR. * & MFR. DESIGNA	ATION	STOCK NO.	#		
R201	Resistor: fixed, composition, 680,000 ohms $\pm 10\%$, 1/2 W	в*	23-680K	1		
R202	Resistor: fixed, composition, 1.2 megohms ±10%, 1/2 W	В*	23-1.2M	2		
R203	Resistor: fixed, composition, 39,000 ohms $\pm 5\%$, 2 W	в*	25-39K-5	3		
R204, 205	Resistor: fixed, composition, 13,000 ohms $\pm 5\%$, 1 W	в*	24-13K-5	2		
R206, 207	Same as R203					
R208	Same as R202					
R209	Part of S202: not separately replaceable					
R210	Resistor: fixed, deposited carbon, 2 megohms ±1%, 1/2 W	NN*	33-2M	1		
R211	Resistor: fixed, composition, 47 ohms $\pm 10\%$, 1/2 W	B*	23-47	1		
R212	Resistor: fixed, composition, $3600 \text{ ohms } \pm 5\%, 1/2 \text{ W}$	в*	23-3600-5	1		
R213	Resistor: fixed, composition, 3300 ohms ±10%, 1/2 W	в*	23-3300	1		
R214	Resistor: fixed, composition, 30,000 ohms $\pm 5\%$, 2 W	в*	25-30K-5	1		
R215	Resistor: fixed, deposited carbon, 216,300 ohms ±1%, 1/2 W	NN*	33-216. 3K	1		
R216	Resistor: fixed, deposited carbon, 252,000 ohms $\pm 1\%$, 1/2 W	NN*	33-252K	2		
R217	Resistor: fixed, deposited carbon, 479,000 ohms ±1%, 1 W	NN*	31-479К	1		
R218	Part of S201: not separately replaceable					
R219	Resistor: fixed, composition, 120,000 ohms $\pm 10\%$, 1/2 W	в*	23-120K	1		
R22 0	Resistor: variable, composition, 100,000 ohms $\pm 30\%$, 1/4 W	BO [‡]	210-138	2		

TABLE OF REPLACEABLE PARTS
R232: Change to resistor, fixed, composition, 1000 ohms <u>+5%</u>, ½ W; -hp- Stock No. 23-1000-5, Mfr., B

CIRCUIT REF.	DESCRIPTION, MFR. * & MFR. DESIGNATION	STOCK NO.	#		
R221	Resistor: fixed, composition, 27,000 ohms ±10%, 1/2 W B*	23-27K	3		
R222	Resistor: fixed, composition, 7.5 megohms ±5%, 1/2 W B*	23-7.5M-5	1		
R223	Resistor: fixed, composition, 16,000 ohms ±5%, 1 W B*	24-16K-5	1		
R224	Resistor: fixed, composition, 62,000 ohms ±5%, 2WB*	25-62K-5	1		
R225, 226	Same as R47				
R227	Same as R216				
R228	Resistor: fixed, deposited carbon, 284,000 ohms ±1%, 1/2 W NN*	33-284K	1		
R229	Resistor: variable, composition, linear taper 500,000 ohms ±30%, 1/4 W BO*	210-146	2		
R230	Resistor: fixed, composition, 390,000 ohms ±10%, 1/2 W B*	23-390K	1		
R231	Resistor: fixed, composition, 270,000 ohms ±10%, 1/2 W B*	23-270K	1		
R232	Resistor: fixed, deposited carbon,680 ohms $\pm 1\%$, $1/2$ WNN*	33-680	1		
R233	Resistor: fixed, composition, 33,000 ohms ±10%, 1 W B*	24-33K	1		
R234	Resistor: fixed, composition, 82,000 ohms ±10%, 2 W B*	25-82K	1		
R235	Resistor: fixed, composition, 470 ohms ±10%, 1/2 W B*	23-470	2		
R236	Same as R221				
R237	Same as R235				
R238	Resistor: fixed, composition, 68,000 ohms ±10%, 2 W B*	25-68K	1		
R239	Same as R8				
R240	Same as R229				

TABLE OF REPLACEABLE PARTS

* See "List of Manufacturers Code Letters For Replaceable Parts Table".

Total quantity used in the instrument.

Model 130B

CIRCUIT REF.	DESCRIPTION, MFR. * & MFR. DESIGNATION	STOCK NO.	#	(act)		
R241	Resistor: fixed, deposited carbon, 2.163 megohms ±1%, 1 W NN*	31-2.163M	1		iner .	
R242	Resistor: fixed, composition, 22 megohms ±10%, 1/2 W B*	23-22M	1		in the	150
R243	Resistor: fixed, composition, 3 megohms ±5%, 1/2 W B*	23-3M-5	1			
R244	Same as R25	1.5		1		
R245	Resistor: fixed, composition, 1200 ohms ±10%, 1/2 W B*	23-1200	1			
R246	Same as R11					
R247 thru R249	These circuit references not assigned	44. 200 s				
R250	Resistor: fixed, deposited carbon, 3.6 megohms ±1%, 1/2 W NN*	33-3.6M	1	10.2		
R251, 252	Resistor: fixed, deposited carbon, 1.2 megohms ±1%, 1/2 W NN*	33-1.2M	1			
R253	Resistor: fixed, deposited carbon, 36 megohms $\pm 1\%$, 2 W NN*	32-36M	1			
R254, 255	Resistor: fixed, deposited carbon, 12 megohms ±1%, 1 W NN*	31-12M	2			
R256	Resistor: fixed, deposited carbon, 1.15 megohms ±1%, 1/2 W NN*	33-1.15M	1		27	
R257, 258	Resistor: fixed, deposited carbon, 376,000 ohms ±1%, 1/2 W NN*	33-376К	2			
R259 thru R261	Resistor: variable, composition, 20,000 ohms ±20%, 1/4 W BO*	210-136	3	-		
R262	Resistor: variable, composition, linear taper, 50,000 ohms ±30%, 1/3 W BO*	210-153	1	al.	10	199
R263	Resistor: fixed, composition, 22,000 ohms ±10%, 1/2 W B*	23-22K	2			area.
R264	This circuit reference not assigned					-

TABLE OF REPLACEABLE PARTS

CIRCUIT REF.	DESCRIPTION, MFR. * & MFR. DESIGNATION	STOCK NO.	#	and the	
R265	Resistor: fixed, composition, 3.9 megohms ±10%, 1/2 W B*	23-3.9M	1	a Shi	
R266 thru R300	These circuit references not assigned	Jul 1			
		C U SP S	120		
R301	Same as R145				
R302	Resistor: fixed, wirewound, 4000 ohms ±5%, 20 W S*	27-7	1		
R303	Same as R220				
R304	Resistor: fixed, composition, 2.2 megohms ±10%, 1 W B*	24-2.2M	2		
R305	Resistor: fixed, composition, 1000 ohms ±10%, 1/2 W B*	23-1000	6	Mind.	
R306	Resistor: fixed, composition, 68,000 ohms ±10%, 1 W B*	24-68K	1	10	
R307	Resistor: fixed, composition, 10,000 ohms ±10%, 1/2 W B*	23-10K	2		
R308	Resistor: fixed, composition, 1 megohm ±10%, 1 W B*	24-1M	2	di ministra	
R309	Same as R305				
R310	Resistor: fixed, deposited carbon, 1.13 megohms ±1%, 1/2 W NN*	33-1.13M	1		
R311	Resistor: fixed, deposited carbon, 1.39 megohms ±1%, 1/2 W NN*	33-1.39M	2		
R312	Resistor: fixed, deposited carbon, 2.84 megohms ±1%, 1 W NN*	31-2.84M	1		
R313	Resistor: fixed, wirewound, 20,000 ohms ±10%, 10 W S*	26-77	1		
R314	Resistor: fixed, wirewound, 1500 ohms ±10%, 20 W S*	27-46	1		
R315	Same as R308				- 1
R316	Same as R305				

TABLE OF REPLACEABLE PARTS

Model 130B

CIRCUIT REF.	DESCRIPTION, MFR. * & MFR. DESIGNATION	STOCK NO.	#	horas		
R317	Resistor: fixed, composition, 22,000 ohms ±10%, 1 W B*	24-22K	1	544A	gert [1.165-1
R318	Same as R221				tenal a	áuc.
R319	Resistor: fixed, composition, 1.2 megohms ±10%, 1 W B*	24-1.2M	2	- / -		1
R320	Same as R305		100	11.9	1.90	122
R321	Resistor: fixed, deposited carbon,1.75 megohms $\pm 1\%$, $1/2$ WNN*	33-1.75M	1		635	
R322	Resistor: fixed, deposited carbon, 2.52 megohms ±1%, 1/2 WNN*	33-2.52M	1	-	an a	
R323	Resistor: fixed, composition, 680,000 ohms ±10%, 1 W B*	24-680K	2	-	ine (
R324	Resistor: fixed, composition, 100,000 ohms ±10%, 1 W B*	24-100K	1	100		10253
R325	Same as R305					100000
R326	Resistor: fixed, composition, 12,000 ohms ±10%, 1 W B*	24-12K	1			
R327	Same as R263			0.0		
R328	Resistor: fixed, composition, 4700 ohms ±10%, 1/2 W B*	23-4700	1	84 1		
R329	Same as R319	1. 22.13		-		
R330	Same as R305		-			
R331	Same as R5	1.1.1.1			l	
R332	Resistor: variable, composition, linear taper, 250,000 ohms ±30%, 1/4 W BO*	210-194	1	1.50		
R333	Same as R311		2.0			1 1
R334	Resistor: variable, composition, linear taper, 1 megohm ±30%, 1/4 W BO*	210-139	1	_		
R335	Resistor: fixed, deposited carbon, 20 megohms ±10%, 1 W AV*	26-84	1		1.5	
						200

TABLE OF REPLACEABLE PARTS

Sect. V Page 18

CIRCUIT REF.	DESCRIPTION, MFR. * & MFR. DESIGNATION	STOCK NO.	#		
R336	Resistor: variable, composition, linear taper, 5 megohms ±30%, 1/2 W I*	210-159	1		
R337	Same as R5		-	Sec.	
R338	Resistor: fixed, composition, 27,000 ohms ±10%, 1/2 W B*	23-27K	1		
R339	Same as R304		1	11.54	
R340	This circuit reference not assigned				
R341	Same as R307	and the second second	1.00		
R342	This circuit reference not assigned	1951 1	10.6	10	
R343	Resistor: variable, composition, 1 megohm ±30%, 0.2 W BO*	210-118	1		-
R344, 345	These circuit references not assigned	Particular.		11	
R346	Resistor: fixed, deposited carbon, 29 megohms ±10%, 1 W AV*	26-85	1		
R347	This circuit reference not assigned			13.11	
R348	Resistor: fixed, composition, 270,000 ohms ±10%, 1 W B*	24-270К	1		
R349	Same as R46				
R350	Resistor: fixed, composition, 47,000 ohms ±10%, 1/2 W B*	23-47K	2		
R351	Resistor: fixed, composition, 1000 ohms ±10%, 1/2 W B*	23-1000	2		
R352	Resistor: fixed, composition, 2000 ohms ±5%, 1/2 W B*	23-2000-5	1	-	
R353	Resistor: fixed, composition, 510 ohms ±5%, 1 W B*	24-510-5	1		
R354	Resistor: fixed, composition, 2700 ohms ±10%, 1/2 W B*	23-2700	1	11260	
R355, 356	Resistor: fixed, composition, 68,000 ohms ±10%, 1/2 W B*	23-68K	2		

TABLE OF REPLACEABLE PARTS

CIRCUIT REF.	DESCRIPTION, MFR. * & MFR. DESIGNATION	Description of the stock no.	#	
R357	Resistor: variable, wirewound, 25 ohms ±10%, 2 W, linear taper BO*	210-188	1	
R358	Same as R323	Sugar the	initia manification to	
R359	Resistor: fixed, deposited carbon, 750,000 ohms ±1%, 1/2 W NN*	33-750K	1	in a
R360	Resistor: fixed, composition, 10 megohms ±10%, 1/2 W B*	23-10M	1 1 1 1 1 - A - A - A - A - A - A - A -	
R361	Resistor: fixed, composition, 56,000 ohms ±10%, 1 W B*	24-56K	1	
R362	Resistor: fixed, composition, 150,000 ohms ±10%, 1 W B*	24-150K	1.000 C	
R363	This circuit reference not assigned	COLUMN COLUMN		
R364	Resistor: fixed, composition, 82,000 ohms ±10%, 1/2 W B*	23-82K	1	
R365	Resistor: variable, composition, linear taper, 50,000 ohms $\pm 20\%$, G*	210-18	1	
R366	Resistor: fixed, composition, 33,000 ohms ±10%, 1/2 W B*	23-33K	1	
R367	This circuit reference not assigned			
R368	Resistor: fixed, composition, 1 megohm ±10%, 1/2 W B*	23-1M	1	
R369	Same as R350	1.4.N-4. VI	Concreanding 1 St.	
R370	Same as R351	SCAT!	are being at	
R371	Resistor: fixed, composition, 2.2 megohms ±10%, 1/2 W B*	23-2.2M	1	
S1	Switch, toggle: SPST Vertical AC-DC Switch D*	310-11	1	
S2	Vertical Sensitivity Switch Assembly HP*	130B-95A	1	
S3	Part of R18: not separately replaceable			N.
1				

TABLE OF REPLACEABLE PARTS

Sect. V Page 20

TABLE	OF	REPLACEABLE	PARTS

CIRCUIT REF.	DESCRIPTION, MFR. * & MFR. DESIGNATION	STOCK NO.	#		1	
				1		
S4 thru S100	These circuit references not assigned	1453				
S101	Switch, rotary: 5 section, 16 pos.	310-236	1	3		
S102	Horizontal Sensitivity Switch Assembly HP*	130B-95B	1		1.2	
S103	Part of R118: not separately replaceable				1	- 1
S104 thru S200	These circuit references not assigned			2001	R	
S201	Sync Switch Assembly HP*	130B-19H	1	100		
S202	Trigger Level Switch Assembly HP*	130B-19G	1			ļ
S203	Part of R218: not separately replaceable					
S204	Sweep Time/CM Switch Assembly HP*	130B-19J	1			
S205	Part of R262: not separately replaceable	Print mark				
S206 thru S300	These circuit references not assigned				1.0	
S301	Part of R357: not separately replaceable				¥.	
		Constant of				ĺ
т1	Transformer, power HP*	910-162	1			
V1, 2	Tube, electron: 12AU7 ZZ*	212-12AU7	4			
V3	Tube, electron: 12AT7 ZZ*	212-12AT7	3			
V4	Tube, electron: 6DJ8 ZZ*	212-6DJ8	3			
V5 thru V100	These circuit references not assigned					
V101, 102	Same as V1					
V103	Same as V3					
V104 XV104	Same as V4					

Model 130B

CIRCUIT REF.	DESCRIPTION	N, MFR. * &MFR. DESK	GNATION	STOCK NO.	#	and the second	(C	
V105 thru V200	These circuit re	ferences not assigned		1999 - C. 1999 - C. 1999		nere e		
V201	Same as V4			57960				
V202	Same as V3			2		1		
V203	Tube, electron:	6U8	ZZ*	212-6U8	1			
V204	Tube, electron:	6C4		212-6C4	1		1	
V205	Tube, electron:	12AL5		212-12AL5	1			
V206	Tube, electron:	6AW8A		212- 6AW8A	1			
V207	Tube, electron:	12AX7	ZZ*	212-12AX7	1	01	- 1	
V208 thru V300	These circuit re	ferences not assigned				deal 1	ind X 1	
V3 01	Tube, electron:	12B4A	ZZ*	212-12B4A	3	100	Dan I.	
V302	Tube, electron:	6AU6	ZZ*	212-6AU6	1	mer	\$9.77	
V303	Same as V301						in the	1411
V304	Tube, electron:	6BH6	ZZ*	212-6BH6	2			10
V305	Tube, electron:	6X4	ZZ*	212-6X4	1	é. an		
V306	Same as V301							
V307	Same as V304							
V308	Tube, electron:	5651	ZZ *	212-5651	1			
V309	Tube, electron, (Final numbe denotes phos	cathode-ray type r in tube designation phor type)	ZZ*	212-5AQP-	1			
V310, 311	Tube, electron:	5642	ZZ*	212-5642	2			
V312	Tube, electron:	12AU7	ZZ*	212-12AU7	1			
V313	Tube, electron:	6AQ5	ZZ*	212-5642	1			
V314	Tube, electron:	6DJ8	ZZ*	212-6DJ8	1			
							iú	

TABLE OF REPLACEABLE PARTS

Sect. V Page 22

....

CIRCUIT REF.	DESCRIPTION, MFR. * & MFR. DESIGNAT	ION	STOCK NO.	#	lint.	
	MISCELLANEOUS	***	1904 904	1		
	CRT bezel	HP*	120A-20A	1		
	Filter, light: amber	HPT	120A-03A	1		
111	Filter, light: blue	HP*	120A-83B	1		
	Filter, light: green	HP*	120A-83G	1		
x 1	Fuseholder	T*	140-16	1		
	Graticule	HP*	120A-83C	1		
	Holder, rectifier	T*	140-91	2		
	Jewel	HP*	145-23A	1		
	Knob: FOCUS, INTENSITY, SCALE LIGHT	HP*	G-74D	3		
	Knob: VERT.POS., HORIZ. POS.	HP*	G-74G	2		
	Knob: TRIGGER LEVEL	HP*	G-74L	1		
	Knob: VERT. SENSITIVITY, HORIZ. SENSITIVITY, SYNC TIME, SWEEP TIME	HP*	G-74Q	4		
	Knob: TRIGGER SLOPE	HP*	G-74AT	1		
	Knob: VERNIER, VERT. SENSITIVITY, HORI SENSITIVITY, SYNC TIME, SWEEP TIME	Z. HP*	G-74AU	4		
	Knob: VERT. and HORIZ. DC BAL	HP*	G-74BJ	2		
	Lampholder Assembly	JJ*	145-27	1		
	Pilot Light Assembly	AD*	145-25	1		
Sec. 1		1. 7	1- m-11			L

TABLE OF REPLACEABLE PARTS

LIST OF CODE LETTERS USED IN TABLE OF REPLACEABLE PARTS TO DESIGNATE THE MANUFACTURERS

CODE LETTER MAN

MANUFACTURER

Aerovox Corp. A В Allen-Bradley Co. С Amperite Co. D Arrow, Hart & Hegeman ε Bussman Manufacturing Co. F Carborundum Co. G Centralah Cinch-Jones Mfg. Co. H Hewlett-Packard Co. HP Clarostat Mfg. Co. Т Cornell Dubilier Elec. Co. J Κ Hi-Q Division of Aerovox Erie Resistor Corp. L Fed. Telephone & Radio Corp. м N General Electric Co. 0 General Electric Supply Corp. Ρ Girard-Hopkins Q Industrial Products Co. R International Resistance Co. Lectrohm Inc. S Ŧ Littlefuse Inc. υ Maguire Industries Inc. ٧ Micamold Radio Corp. W Oak Manufacturing Co. P. R. Mallory Co., Inc. Х Y Radio Corp. of America Sangamo Electric Co. Ζ AA Sarkes Tarzian Signal Indicator Co. BB CC Sprague Electric Co. DD Stackpole Carbon Co. EE Sylvania Electric Products Co. FF Western Electric Co. GG Wilkor Products, Inc. HH Amphenol 11 Dial Light Co. of America JJ Leecraft Manufacturing Co. KK Switchcraft, Inc. LL Gremar Manufacturing Co. Carad Corp. MM NN Electra Manufacturing Co. 00 Acro Manufacturing Co. PP Alliance Manufacturing Co. QQ Arco Electronics, Inc. RR Astron Corp. SS Axel Brothers Inc. TT Belden Manufacturing Co. UU Bird Electronics Corp. ٧V Barber Colman Co. ww Bud Radio Inc. Allen D. Cardwell Mfg. Co. XX YY Cinema Engineering Co. Any brand tube meeting ZZ **RETMA** standards. AB Corning Glass Works AC Dale Products, Inc. AD The Drake Mfg. Co. AE Elco Corp. AF Hugh H. Eby Co. AG Thomas A. Edison, Inc. Fansteel Metallurgical Corp. AH General Ceramics & Steatite Corp. AL

The Gudeman Co.

AJ

ADDRESS

New Bedford, Mass. Milwaukee 4, Wis. New York, N.Y. Hartford, Conn. St. Louis, Mo. Niagara Falls, N.Y. Milwaukee I, Wis. Chicago 24, Ill. Palo Alto, Calif. Dover, N. H. South Plainfield, N. J. Olean, N.Y. Erie 6. Pa. Clifton, N. J. Schenectady 5, N.Y. San Francisco, Calif. Oakland, Calif. Danbury, Conn. Philadelphia 8, Pa. Chicago 20, Ill. Des Plaines, Ill. Greenwich, Conn. Brooklyn 37, N.Y. Chicago 10, Ill. Indianapolis, Ind. Harrison, N. J. Marion, Ill. Bloomington, Ind. Brooklyn 37, N.Y. North Adams, Mass. St. Marys, Pa. Warren, Pa. New York 5, N.Y. Cleveland, Ohio Chicago 50, Ill. Brooklyn 37, N.Y. New York, N.Y. Chicago 22, Ill. Wakefield, Mass, Redwood City, Calif. Kansas City, Mo. Columbus 16, Ohio Alliance, Ohio New York 13, N.Y. East Newark, N. J. Long Island City, N.Y. Chicago 44, Ill. Cleveland 14, Ohio Rockford, Ill. Cleveland 3. Ohio Plainville, Conn. Burbank, Calif.

Corning, N. Y. Columbus, Neb. Chicago 22, III. Philadelphia 24, Pa. Philadelphia 44, Pa. West Orange, N. J. North Chicago, III. Keasbey, N. J. Sunnyvale, Calif.

cw

Good-All Electric Mfg. Co.

CODE LETTER MANUFACTURER

AK Hammerlund Mfg. Co., Inc. Industrial Condenser Corp. AL AM Insuline Corp. of America Jennings Radio Mfg. Corp. AN E. F. Johnson Co. AO Lenz Electric Mfg. Co. AP AQ Micro-Switch AR Mechanical Industries Prod. Co. AS Model Eng. & Mfg., Inc. The Muter Co. AT AU Ohmite Mfg. Co. Resistance Products Co. AV AW Radio Condenser Co. AX Shallcross Manufacturing Co. AY Solar Manufacturing Co. AZ Sealectro Corp. BA Spencer Thermostat BC Stevens Manufacturing Co. BD Torrington Manufacturing Co. BE Vector Electronic Co. BF Weston Electrical Inst. Corp. BG Advance Electric & Relay Co. BH E. I. DuPont BI Electronics Tube Corp. BJ Aircraft Radio Corp. BK Allied Control Co., Inc. BL Augat Brothers, Inc. BM **Carter Radio Division** BN **CBS Hytron Radio & Electric** BO Chicago Telephone Supply Henry L. Crowley Co., Inc. BP Curtiss-Wright Corp. BQ BR Allen B. DuMont Labs BS Excel Transformer Co. BT General Radio Co. BU Hughes Aircraft Co. RV International Rectifier Corp. BW James Knights Co. ΒX Mueller Electric Co. Precision Thermometer & Inst. Co. BY ΒZ Radio Essentials Inc. CA Raytheon Manufacturing Co. CB Tung-Sol Lamp Works, Inc. CD Varian Associates CE Victory Engineering Corp. CF Weckesser Co. CG Wilco Corporation СН Winchester Electronics, Inc. CI Malco Tool & Die CJ Oxford Electric Corp. CK Camloc-Fastener Corp. CL George K. Garrett СМ Union Switch & Signal Radio Receptor CN co Automatic & Precision Mfg. Co. CP Bassick Co. CQ Birnbach Radio Co. CR **Fischer Specialties** Telefunken (c/o MVM, Inc.) CS СТ Potter-Brumfield Co. Cυ Cannon Electric Co. C٧ Dynac, Inc.

ADDRESS

New York I, N.Y. Chicago 18, Ill. Manchester, N. H. San Jose, Calif. Waseca, Minn. Chicago 47, Ill. Freeport, Ill. Akron 8, Ohio Huntington, Ind. Chicago 5, Ill. Skokie, Ill. Harrisburg, Pa. Camden 3, N. J. Collingdale, Pa. Los Angeles 58, Calif. New Rochelle, N.Y. Attleboro, Mass. Mansfield, Ohio Van Nuys, Calif. Los Angeles 65, Calif. Newark 5, N. J. Burbank, Calif. San Francisco, Calif. Philadelphia 18, Pa. Boonton, N. J. New York 21, N.Y. Attleboro, Mass. Chicago, III. Danvers, Mass. Elkhart, Ind. West Orange, N. J. Carlstadt, N. J. Clifton, N. J. Oakland, Calif. Cambridge 39, Mass. Culver City, Calif. El Segundo, Calif. Sandwich, III. Cleveland, Ohio Philadelphia 30, Pa. Mt. Vernon, N.Y. Newton, Mass. Newark 4, N. J. Palo Alto, Calif. Union, N. J. Chicago 30, Ill. Indianapolis, Ind. Santa Monica, Calif. Los Angeles 42, Calif. Chicago 15, Ill. Paramus, N. J. Philadelphia 34, Pa. Swissvale, Pa. New York 11, N.Y. Yonkers, N.Y. Bridgeport 2, Conn. New York 13, N.Y. Cincinnati 6, Ohio New York, N.Y. Princeton, Ind. Los Angeles, Calif. Palo Alto, Calif. Ogallala, Nebr.

CLAIM FOR DAMAGE IN SHIPMENT

The instrument should be tested as soon as it is received. If it fails to operate properly, or is damaged in any way, a claim should be filed with the carrier. A full report of the damage should be obtained by the claim agent, and this report should be forwarded to us. We will then advise you of the disposition to be made of the equipment and arrange for repair or replacement. Include model number and serial number when referring to this instrument for any reason.

WARRANTY

Hewlett-Packard Company warrants each instrument manufactured by them to be free from defects in material and workmanship. Our liability under this warranty is limited to servicing or adjusting any instrument returned to the factory for that purpose and to replace any defective parts thereof. Klystron tubes as well as other electron tubes, fuses and batteries are specifically excluded from any liability. This warranty is effective for one year after delivery to the original purchaser when the instrument is returned, transportation charges prepaid by the original purchaser, and when upon our examination it is disclosed to our satisfaction to be defective. If the fault has been caused by misuse or abnormal conditions of operation, repairs will be billed at cost. In this case, an estimate will be submitted before the work is started.

If any fault develops, the following steps should be taken:

1. Notify us, giving full details of the difficulty, and include the model number and serial number. On receipt of this information, we will give you service data or shipping instructions.

2. On receipt of shipping instructions, forward the instrument prepaid, to the factory or to the authorized repair station indicated on the instructions. If requested, an estimate of the charges will be made before the work begins provided the instrument is not covered by the warranty.

SHIPPING

All shipments of Hewlett-Packard instruments should be made via Truck or Railway Express. The instruments should be packed in a strong exterior container and surrounded by two or three inches of excelsior or similar shock-absorbing material.

DO NOT HESITATE TO CALL ON US

HEWLETT-PACKARD COMPANY Laboratory Instruments for Speed and Accuracy 275 PAGE MILL ROAD CABLE PALO ALTO. CALIF. U.S.A. "HEWPACK"

